

Datasheet: VMA00059KT

ZAP 70 ANTIBODY WITH CONTROL LYSATE
ZAP 70
Purified
PrecisionAb Monoclonal
lgG1
2 Westerns

### **Product Details**

### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting				1/1000

**PrecisionAb antibodies have been extensively** <u>validated for the western blot</u> <u>application.</u> The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	20µl Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN <sub>3</sub> )
Immunogen	Human ZAP 70 amino acids residues 280-309
External Database Links	UniProt:

P43403 Related reagents

### **Entrez Gene:**

7535 ZAP70 Related reagents

Synonyms	SRK				
Specificity	<b>Mouse anti Human ZAP 70 antibody</b> recognizes human ZAP-70, a 70 kDa tyrosine protein kinase which associates with the T-cell receptor (TCR) zeta chain and undergoes phosphorylation following TCR stimulation.				
	ZAP-70 is primarily expressed in T-cells and natural killer cells, where it plays a key role in T-cell receptor (TCR) signalling. Recent studies suggest that ZAP-70 is also expressed in a population of normal immature B-cells and in B-cells from a sub-set of patients with chronic lymphocytic leukaemia.				
	Mouse anti Human ZAP 70 antibody recognizes ZAP 70 as a 70 kDa band in western blotting on Jurkat and MOLT-4 cell lysates.				
Western Blotting	Anti ZAP 70 antibody detects a band of approximately 70 kDa in Jurkat cell lysates.				
Instructions For Use	Please refer to the <u>PrecisionAb western blotting protocol.</u> For additional information on secondary antibody dilution and exposure time see product web page.				
Lysate Composition	400μg Jurkat lysate lyophilized in RIPA buffer.				
Lysate Reconstitution	- If using DDT reconstitute the lyophilized lysate with 190μl DI H <sub>2</sub> O, add 200μl 2x Laemmli Sample Buffer and 10μl 2M DTT If using BME reconstitute the lyophilized lysate with 180μl DI H <sub>2</sub> O, add 200μl 2x Laemml Sample Buffer and 20μl BME. Heat at 95°C for 5 minutes. For 10 well mini gels load 25μl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.				
Storage	Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.				
	Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.				
Guarantee	As supplied, 12 months from date of despatch.				
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.				
Health And Safety Information	Material Safety Datasheet documentation #10077 #10561 available at: Antibody (10077): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10077.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10077.pdf</a> Lysate Material (10561): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf</a>				

# **Related Products**

### **Recommended Secondary Antibodies**

Goat Anti Mouse IgG (H/L) (STAR207...) HRP

## **Recommended Negative Controls**

#### MOUSE IgG1 NEGATIVE CONTROL (MCA928)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

From March 15, 2021, we will no longer supply printed datasheets with our products. Look out for updates on how to access your digital version at bio-rad-antibodies.com 'M351729:190318'

Printed on 09 Feb 2021

© 2021 Bio-Rad Laboratories Inc | Legal | Imprint