

Datasheet: 0100-0242

**BATCH NUMBER 159578**

<b>Description:</b>	MOUSE ANTI HUMAN CD121a
<b>Specificity:</b>	CD121a
<b>Other names:</b>	INTERLEUKIN 1 RECEPTOR TYPE 1
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	35730 (49/20)
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	<5% Trehalose
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Recombinant peptide corresponding to the extracellular domain of human CD121a

expressed in insect cell line Sf21.

---

**External Database****Links****UniProt:**[P14778](#)[Related reagents](#)**Entrez Gene:**[3554](#)

IL1R1

[Related reagents](#)

---

**Synonyms**

IL1R, IL1RA, IL1RT1

---

**RRID**

AB\_617210

---

**Specificity**

**Mouse anti human CD121a, clone 35730**, recognizes interleukin-1 receptor type 1 (IL1R1). IL1RA/CD121a is widely expressed in a number of cell types including T cells and fibroblasts ([Sims \*et al.\*](#)) and is involved in the activation of NF-kappa-B ([Slack \*et al.\*](#)).

---

**Further Reading**

1. Sims, J. *et al.* (1989) Cloning the interleukin 1 receptor from human T cells. [Proc Natl Acad Sci U S A. 86: 8946-50.](#)
2. Slack, J. *et al.* (2000) Identification of two major sites in the type I interleukin-1 receptor cytoplasmic region responsible for coupling to pro-inflammatory signaling pathways. [J Biol Chem. 275:4670-8.](#)

---

**Storage**

Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10121 available at:  
<https://www.bio-rad-antibodies.com/SDS/0100-0242>  
10121

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)Goat Anti Mouse IgG (STAR76...) [RPE](#)Goat Anti Mouse IgG (STAR70...) [FITC](#)Rabbit Anti Mouse IgG (STAR13...) [HRP](#)Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

Goat Anti Mouse IgG (STAR77...)

[HRP](#)

Goat Anti Mouse IgG (H/L) (STAR117...)

[Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

**North & South  
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M379777:210413'

**Printed on 15 Apr 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)