

Datasheet: 0151-9104

**BATCH NUMBER 180501**

|                      |                             |
|----------------------|-----------------------------|
| <b>Description:</b>  | GOAT ANTI ADENOVIRUS:Biotin |
| <b>Specificity:</b>  | ADENOVIRUS                  |
| <b>Format:</b>       | Biotin                      |
| <b>Product Type:</b> | Polyclonal Antibody         |
| <b>Isotype:</b>      | Polyclonal IgG              |
| <b>Quantity:</b>     | 1 ml                        |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                    | Yes | No | Not Determined | Suggested Dilution |
|--------------------|-----|----|----------------|--------------------|
| ELISA              | ▪   |    |                | 1/5,000            |
| Western Blotting   | ▪   |    |                | 1/1,000            |
| Immunofluorescence | ▪   |    |                | 1/8,000            |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

|                                       |   |
|---------------------------------------|---|
| <b>Target Species</b>                 | Viral   |
| <b>Product Form</b>                   | Purified IgG conjugated to Biotin - liquid  |
| <b>Buffer Solution</b>                | Phosphate buffered saline   |
| <b>Preservative Stabilisers</b>       | 0.09% Sodium Azide (NaN <sub>3</sub> )  |
| <b>Approx. Protein Concentrations</b> | 1.0 mg/ml   |
| <b>Immunogen</b>                      | Gradient purified, disrupted native adenovirus type 5 virions.                          |
| <b>RRID</b>                           | AB_616736   |
| <b>Specificity</b>                    | <b>Goat anti Adenovirus antibody</b> recognizes numerous Adeno proteins including hexon |

protein, and reacts with Adenovirus types 1, 2, 3, 5, 6, 7a, 8, 31, 40 & 41. Broad reactivity with other Adenovirus types is expected.

Goat anti Adenovirus antibody does not cross-react with Parainfluenza types 1-3, Influenza A & B, RSV, hMPV and Rhinovirus type 16. Goat anti Adenovirus antibody does not react with uninfected Hep-2, LLCMK2 and Vero cells.

Adenoviruses are 80-110 nm diameter, non-enveloped icosahedral viruses that can cause respiratory illness, infantile gastroenteritis and other mild disorders such as conjunctivitis. Features of Adenovirus have made it a useful tool for gene transfection and expression in mammalian cells.

---

**References**

1. Mallery, D.L. *et al.* (2010) Antibodies mediate intracellular immunity through tripartite motif-containing 21 (TRIM21). [Proc Natl Acad Sci U S A. 107 \(46\): 19985-90.](#)
2. Chappell, C.L. *et al.* (2016) Adenovirus 36 Antibody Detection: Improving the Standard Serum Neutralization Assay. [J Virol Methods. Nov 5. pii: S0166-0934\(16\)30405-0. \[Epub ahead of print\]](#)

---

**Storage**

Store at +4°C or at -20°C if preferred.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/0151-9104>  
10040

---

**Regulatory**

For research purposes only

---

**North & South America**

Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M362759:200528'

**Printed on 01 Mar 2024**