

Datasheet: 0178-0149

**BATCH NUMBER 169201**

<b>Description:</b>	MOUSE ANTI HUMAN ACTH
<b>Specificity:</b>	ACTH
<b>Other names:</b>	ADRENOCORTICOTROPIC HORMONE
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	57 (BGN/1388/66)
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1/100 - 1/200
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1) This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase.**

### Target Species

Human

### Species Cross Reactivity

Reacts with: Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	Synthetic peptide corresponding to the amino acids 1 - 24 (SYSMEHFRWGKPVGKKRRPVKVYP) of the human ACTH conjugated to Keyhole Limpet Haemocyanin.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01189</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">5443</a>    POMC    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_2166032
<b>Specificity</b>	<p><b>Mouse anti Human ACTH antibody (Clone 57, BGN/1388/66)</b> recognizes 1-24 ACTH (Synacthen) which, as a synthetic analogue of naturally-occurring Adrenocorticotrophic Hormone, can be used to measure adrenal reserve and for the diagnosis of adrenal insufficiency by acute adrenocortical stimulation.</p> <p>ACTH, released from the anterior pituitary gland in response to corticotropin-releasing hormone from the hypothalamus, acts on the adrenal cortex to stimulate the production of corticosteroids such as cortisol, involved in the response to stress. Administration of 1-24 ACTH increases blood pressure, believed to be attributed to an increase in ACTH-stimulated cortisol secretion, in association with increased cardiac output.</p> <p><b>Mouse anti Human ACTH antibody</b> does not recognize 17-39 ACTH (CLIP). It does cross-reacts with rat N-terminal ACTH.</p>
<b>References</b>	1. Zhou, M <i>et al.</i> (2016) Localization of ATP-sensitive K <sup>+</sup> channel subunits in rat pituitary gland. <a href="#">Arch Histol Cytol. 76 (1): 53-66.</a>
<b>Further Reading</b>	1. Hatakeyama, H. <i>et al.</i> (2000) Functional adrenocorticotrophic hormone receptor in cultured human vascular endothelial cells : possible role in control of blood pressure. <a href="#">Hypertension. 36 (5): 862-5.</a>

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/0178-0149">https://www.bio-rad-antibodies.com/SDS/0178-0149</a> 10040
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">Alk. Phos.</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M405238:220916'

Printed on 29 Aug 2024