

Datasheet: 0300-0135

BATCH NUMBER 156048

Description:	GOAT F(ab) ANTI RABBIT IgG (H/L)
Specificity:	IgG (H/L)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Rabbit
Product Form	F(ab) fragment of IgG - liquid
Preparation	Purified F(ab) prepared from monospecific affinity-purified IgG by papain digestion followed by chromatographic separation.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	F(ab) concentration 1.0mg/ml
Immunogen	Rabbit IgG whole molecule.
External Database Links	UniProt:

P01840	Related reagents
P01870	Related reagents
P01839	Related reagents
P01841	Related reagents
P01838	Related reagents
P01847	Related reagents

Entrez Gene:

100009097	LOC100009097	Related reagents
100422811	IGKC2	Related reagents

RRID	AB_618738
Specificity	Goat F(ab') anti Rabbit IgG (H+L) antibody recognizes the heavy and light chains of rabbit IgG. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-goat serum. No reaction was observed against anti-papain or anti-goat IgG Fc.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/0300-0135
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (H/L) (602001...) [Alk. Phos.](#), [FITC](#)
 Donkey Anti Goat IgG (H/L) (642001...) [FITC](#), [HRP](#), [RPE](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
 'M362805:200528'

Printed on 20 Aug 2025