

Datasheet: 0300-0218

BATCH NUMBER 171038

Description:	RABBIT ANTI STREPTOCOCCUS PNEUMONIAE
Specificity:	STREPTOCOCCUS PNEUMONIAE
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
ELISA			▪	
Western Blotting			▪	
Immunofluorescence	▪			1/500

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 5.0 mg/ml
Immunogen	Whole cell blend of numerous serotypes.

RRID AB_620715

Specificity Rabbit anti *Streptococcus pneumoniae* antibody recognizes *Streptococcus pneumoniae* and reacts with types 3, 4, 6, 7, 9, 14, 18, 19 and 23.

S. pneumoniae is a non-motile gram-positive bacteria that is found in the human upper respiratory tract and often as diplococci (pairs of cocci). Rabbit anti *Streptococcus pneumoniae* antibody is not absorbed and may react with related microorganisms.

References 1. Doyle, C.R. & Pirofski, L.A. (2016) Reduction of *Streptococcus pneumoniae* Colonization and Dissemination by a Nonopsonic Capsular Polysaccharide Antibody. *MBio.* 7 (1): e02260-15.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/0300-0218>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M419986:230706'

Printed on 29 Jan 2026
