

Datasheet: 0300-1706

**BATCH NUMBER 168983**

<b>Description:</b>	NATIVE CALF ALKALINE PHOSPHATASE
<b>Name:</b>	ALKALINE PHOSPHATASE
<b>Other names:</b>	AP
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**Product Form** Purified AP from calf intestine - liquid

**Preparation** Crystalline suspension. Free of endonuclease, exonuclease and RNase activities.

**Buffer Solution** 70% Ammonium sulphate.

**Preservative Stabilisers** None present

**Approx. Protein Concentrations** Total protein concentration 15.0 mg/ml

### External Database Links

**UniProt:**

[P19111](#) [Related reagents](#)

**Entrez Gene:**

[280993](#) ALPI [Related reagents](#)

<b>Protein Molecular Weight</b>	110-150 kDa
<b>Activity</b>	4000 U/mg protein. 1 unit liberates 1 µmole of p-nitrophenol per minute at 37°C and pH 9.8 in DEA buffer. Activity is determined by measuring the increase in the optical density at 405nm resulting from the hydrolysis of p-nitrophenyl phosphate.
<b>Instructions For Use</b>	The product may be diluted in sterile distilled water or in dilute buffer
<b>References</b>	1. Winkelman, J. <i>et al.</i> (1972) The clinical usefulness of alkaline phosphatase isoenzyme determinations. <a href="#">Am J Clin Pathol. 57 (5): 625-34.</a>
<b>Storage</b>	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10315 available at: <a href="https://www.bio-rad-antibodies.com/SDS/0300-1706">https://www.bio-rad-antibodies.com/SDS/0300-1706</a> 10315
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M362825:200528'

Printed on 29 Aug 2024