

Datasheet: 0650-0311

BATCH NUMBER 161748

Description:	NATIVE HUMAN APOLIPOPROTEIN A1
Name:	APOLIPOPROTEIN A1
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified apolipoprotein A-1 from human plasma - lyophilized
Reconstitution	Reconstitute with 1.0 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.
Preparation	Purified human Apo A1 prepared by ultracentrifugation, delipidation and gel filtration
Buffer Solution	0.05 M Sodium chloride, 0.01 M Sodium carbonate
Preservative Stabilisers	None present
External Database Links	UniProt: P02647 Related reagents

Entrez Gene:

[335](#) APOA1 [Related reagents](#)

Product Information **Native Human apolipoprotein A1** can be used for coating microplates and as a ligand for immunosorbent preparation.

Apolipoproteins are lipid-binding proteins which enable the transport of dietary lipids for storage, metabolism and secretion. Apolipoprotein A-1 (also known as Apo-A1) plays an important part in the removal of cholesterol from cells.

Protein Molecular Weight 28 kDa

Purity >90% by SDS PAGE

References

1. Tall, A.R. & Small, D.M. (1980) Body cholesterol removal: role of plasma high-density lipoproteins. [Adv Lipid Res. 17: 1-51.](#)
2. Fielding, C.J. (1972) A protein cofactor of lecithin:cholesterol acyltransferase. [Biochem. Biophys. Res. Commun. 46: 1493-1498.](#)

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10449 available at:
<https://www.bio-rad-antibodies.com/SDS/0650-0311>
10449

Donor material tested and found negative for HBsAg, HIV1, HIV2 and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350056:190307'

Printed on 20 Mar 2024