

Datasheet: 101001

BATCH NUMBER 167798

Description:	GOAT ANTI MOUSE IgG/A/M (H/L)
Specificity:	IgG IgA IgM (H/L)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoblotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to Mouse IgG, IgA and IgM were raised by repeated immunisation of goats with pooled mouse IgG, IgA and IgM paraproteins. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present.
Approx. Protein Concentrations	IgG concentration 1mg/ml

External Database

Links

UniProt:

P01864	Related reagents
P01872	Related reagents
P03987	Related reagents
P01873	Related reagents
P01868	Related reagents
P01865	Related reagents
P01869	Related reagents
P01878	Related reagents
P01867	Related reagents
P01863	Related reagents
P01844	Related reagents
P01843	Related reagents
P01845	Related reagents
P01834	Related reagents

Entrez Gene:

380793	Igh-1a	Related reagents
16019	Ighm	Related reagents
16017	Ighg1	Related reagents
16019	Ighm	Related reagents
380793	Igh-1a	Related reagents
380795	AI324046	Related reagents
3514	IGKC	Related reagents
16016	Ighg2b	Related reagents
16017	Ighg1	Related reagents
16061	Igh-VJ558	Related reagents
110786	IgIc2	Related reagents
110787	IgIc3	Related reagents
380793	Igh-1a	Related reagents
433053	LOC433053	Related reagents

Synonyms

Igh-4

RRID

AB_609680

Specificity

Goat anti Mouse IgG/A/M reacts with the heavy and light chains of all major classes of mouse immunoglobulin.

The antibody has been adsorbed against human serum to minimise cross-reactivity with human immunoglobulins.

References

1. Ding, G. *et al.* (2010) Effect of cryopreservation on biological and immunological properties of stem cells from apical papilla. [J Cell Physiol. 22: 415-22.](#)
2. Shao, L. *et al.* (2015) Comparative *In Vitro* and *In Vivo* Studies of Porcine Rotavirus G9P[13] and Human Rotavirus Wa G1P[8]. [J Virol. 90 \(1\): 142-51.](#)
3. Iver, S. *et al.* (2019) Probing BAK and BAX Activation and Pore Assembly with Cytochrome c Release, Limited Proteolysis, and Oxidant-Induced Linkage. [Methods Mol Biol. 1877: 201-16.](#)

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/101001>
10123

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America

Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M420274:230706'

Printed on 01 Mar 2024