

Datasheet: 101005

BATCH NUMBER 155482

Description:	GOAT ANTI MOUSE IgG/A/M (H/L):HRP
Specificity:	IgG IgA IgM (H/L)
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Immunoprecipitation			▪	
Western Blotting	▪			1/2000 - 1/10,000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Product Form

Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation

Antisera to mouse IgG, A, M were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

50% Glycerol

Immunogen

Purified mouse IgG, IgA and IgM.

**External Database
Links**

UniProt:

P01864	Related reagents
P01872	Related reagents
P03987	Related reagents
P01873	Related reagents
P01868	Related reagents
P01865	Related reagents
P01869	Related reagents
P01878	Related reagents
P01867	Related reagents
P01863	Related reagents
P01844	Related reagents
P01843	Related reagents
P01845	Related reagents
P01834	Related reagents

Entrez Gene:

380793	Igh-1a	Related reagents
16019	Ighm	Related reagents
16017	Ighg1	Related reagents
16019	Ighm	Related reagents
380793	Igh-1a	Related reagents
380795	AI324046	Related reagents
3514	IGKC	Related reagents
16016	Ighg2b	Related reagents
16017	Ighg1	Related reagents
16061	Igh-VJ558	Related reagents
110786	Iglc2	Related reagents
110787	Iglc3	Related reagents
380793	Igh-1a	Related reagents
433053	LOC433053	Related reagents

Synonyms Igh-4

RRID AB_609673

Specificity **Goat anti Mouse IgG/A/M** reacts with the heavy and light chains of all major classes of mouse immunoglobulin.

The antibody has been adsorbed against human serum to minimise cross-reactivity with human immunoglobulins.

References 1. Ding, G. *et al.* (2010) Effect of cryopreservation on biological and immunological properties of stem cells from apical papilla. [J Cell Physiol. 22: 415-22.](#)

2. Shao, L. *et al.* (2015) Comparative *In Vitro* and *In Vivo* Studies of Porcine Rotavirus G9P[13] and Human Rotavirus Wa G1P[8]. [J Virol. 90 \(1\): 142-51.](#)
3. Iver, S. *et al.* (2019) Probing BAK and BAX Activation and Pore Assembly with Cytochrome c Release, Limited Proteolysis, and Oxidant-Induced Linkage. [Methods Mol Biol. 1877: 201-16.](#)

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10106 available at: https://www.bio-rad-antibodies.com/SDS/101005 10106
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350058:190307'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)