

Datasheet: 1012701

BATCH NUMBER 155048

Description:	RAT ANTI HUMAN G-CSF
Specificity:	G-CSF
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	BVD13-3A5
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			1ug/ml - 5ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant human granulocyte-colony stimulating factor expressed in <i>E. coli</i> .
External Database Links	UniProt: P09919 Related reagents

Entrez Gene:[1440](#) CSF3 [Related reagents](#)

Synonyms	C17orf33, GCSF
RRID	AB_620005
Specificity	<p>Rat anti Human G-CSF antibody, clone BVD13-3A5 recognizes human granulocyte colony stimulating factor (G-CSF), also known as Pluripoietin. In addition Filgrastim and Lenograstim are used to describe recombinant human G-CSF preparations. Human G-CSF is a 207 amino acid, ~20 kDa O-linked glycoprotein belonging to the IL-6 superfamily. There are potentially three isoforms formed by alternative splicing, Rat anti Human G-CSF antibody, clone BVD13-3A5 binds the canonical isoform-1, binding to other isoforms has not been assessed (UniProt: P09919).</p> <p>G-CSF is produced by activated macrophages, monocytes and neutrophils. In addition it can be produced by fibroblasts and endothelial cells (Zsebo et al. 1988). G-CSF stimulates production of cells belonging to the granulocytic lineage in the bone marrow and is an activating agent for mature neutrophilic granulocytes (Scapini et al. 2003) to produce other cytokines including CD257 (BAFF), operating through binding to the G-CSF receptor CD114.</p> <p>Rat anti Human G-CSF antibody, clone BVD13-3A5 has been successfully evaluated as a capture reagent in a highly sensitive sandwich ELISA assay with biotinylated Rat anti Human G-CSF, clone BVD11-37G10 (1012808) as a detection reagent and recombinant human G-CSF (PHP082B) as antigen.</p>
ELISA	This product may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with 1012808 as the detection reagent.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10123 available at: https://www.bio-rad-antibodies.com/SDS/1012701 10123
Regulatory	For research purposes only

Related Products

[Recommended Secondary Antibodies](#)

Rabbit Anti Rat IgG (STAR16...)	DyLight@800
Rabbit Anti Rat IgG (STAR17...)	FITC
Goat Anti Rat IgG (STAR72...)	HRP
Goat Anti Rat IgG (STAR69...)	FITC
Goat Anti Rat IgG (STAR73...)	RPE
Rabbit Anti Rat IgG (STAR21...)	HRP
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	DyLight@550 , DyLight@650 , DyLight@800
Goat Anti Rat IgG (STAR131...)	Alk. Phos. , Biotin

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M350066:190307'

Printed on 18 Jan 2024