

Datasheet: 1012701

BATCH NUMBER 174208

Description:	RAT ANTI HUMAN G-CSF
Specificity:	G-CSF
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	BVD13-3A5
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			1ug/ml - 5ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Antibody purified from ascites
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant human granulocyte-colony stimulating factor expressed in <i>E. coil</i> .

**External Database
Links**

UniProt:

[P09919](#)

[Related reagents](#)

Entrez Gene:

[1440](#)

CSF3

[Related reagents](#)

Synonyms

C17orf33, GCSF

RRID

AB_620005

Specificity

Rat anti Human G-CSF antibody, clone BVD13-3A5 recognizes human granulocyte colony stimulating factor (G-CSF), also known as Pluripoietin. In addition Filgrastim and Lenograstim are used to describe recombinant human G-CSF preparations. Human G-CSF is a 207 amino acid, ~20 kDa O-linked glycoprotein belonging to the IL-6 superfamily. There are potentially three isoforms formed by alternative splicing, Rat anti Human G-CSF antibody, clone BVD13-3A5 binds the canonical isoform-1, binding to other isoforms has not been assessed ([UniProt: P09919](#)).

G-CSF is produced by activated macrophages, monocytes and neutrophils. In addition it can be produced by fibroblasts and endothelial cells ([Zsebo *et al.* 1988](#)). G-CSF stimulates production of cells belonging to the granulocytic lineage in the bone marrow and is an activating agent for mature neutrophilic granulocytes ([Scapini *et al.* 2003](#)) to produce other cytokines including CD257 (BAFF), operating through binding to the G-CSF receptor CD114.

Rat anti Human G-CSF antibody, clone BVD13-3A5 has been successfully evaluated as a capture reagent in a highly sensitive sandwich ELISA assay with biotinylated Rat anti Human G-CSF, clone BVD11-37G10 ([1012808](#)) as a detection reagent and recombinant human G-CSF as antigen.

ELISA

This product may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with [1012808](#) as the detection reagent.

Storage

Store at -20°C only.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

**Health And Safety
Information**

Material Safety Datasheet documentation #10123 available at:
<https://www.bio-rad-antibodies.com/SDS/1012701>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...)	DyLight®800
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	DyLight®550 , DyLight®650 , DyLight®800
Goat Anti Rat IgG (STAR73...)	RPE
Goat Anti Rat IgG (STAR131...)	Alk. Phos. , Biotin
Rabbit Anti Rat IgG (STAR21...)	HRP
Rabbit Anti Rat IgG (STAR17...)	FITC
Goat Anti Rat IgG (STAR69...)	FITC
Goat Anti Rat IgG (STAR72...)	HRP

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M405695:220916'

Printed on 19 Aug 2026