

Datasheet: 1012801

**BATCH NUMBER 163118**

<b>Description:</b>	RAT ANTI HUMAN G-CSF
<b>Specificity:</b>	G-CSF
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	BVD11-37G10
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			1.0 - 5.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Buffer Solution</b>	Borate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5mg/ml

### External Database Links

**UniProt:**  
[P09919](#)   [Related reagents](#)

**Entrez Gene:**

---

<b>Synonyms</b>	C17orf33, GCSF
-----------------	----------------

---

<b>RRID</b>	AB_618524
-------------	-----------

---

<b>Specificity</b>	<p><b>Rat anti Human G-CSF antibody, clone BVD11-37G10</b> recognizes human granulocyte colony stimulating factor (G-CSF), also known as Pluripoietin. In addition Filgrastim and Lenograstim are used to describe recombinant human G-CSF preparations. Human G-CSF is a 207 amino acid, ~20 kDa O-linked glycoprotein belonging to the IL-6 superfamily. There are potentially three isoforms formed by alternative splicing, Rat anti Human G-CSF antibody, clone BVD13-3A5 binds the canonical isoform-1, binding to other isoforms has not been assessed (<a href="#">UniProt: P09919</a>).</p> <p>G-CSF is produced by activated macrophages, monocytes and neutrophils. In addition it can be produced by fibroblasts and endothelial cells (<a href="#">Zsebo et al. 1988</a>). G-CSF stimulates production of cells belonging to the granulocytic lineage in the bone marrow and is an activating agent for mature neutrophilic granulocytes (<a href="#">Scapini et al. 2003</a>) to produce other cytokines including <a href="#">CD257</a> (BAFF), operating through binding to the G-CSF receptor, CD114.</p> <p>Biotinylated Rat anti Human G-CSF antibody, clone BVD11-37G10 has been successfully evaluated as a detection reagent in a highly sensitive sandwich ELISA assay with biotinylated Rat anti Human G-CSF, clone BVD11-3A5 (<a href="#">1012701</a>) as a capture reagent and recombinant human G-CSF (<a href="#">PHP082B</a>) as antigen.</p>
--------------------	--

---

<b>Storage</b>	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	---

---

<b>Guarantee</b>	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10123 available at: <a href="https://www.bio-rad-antibodies.com/SDS/1012801">https://www.bio-rad-antibodies.com/SDS/1012801</a> 10123
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...)	<a href="#">DyLight®800</a>
Rabbit Anti Rat IgG (STAR17...)	<a href="#">FITC</a>
Goat Anti Rat IgG (STAR72...)	<a href="#">HRP</a>

Goat Anti Rat IgG (STAR69...) [FITC](#)  
Goat Anti Rat IgG (STAR73...) [RPE](#)  
Rabbit Anti Rat IgG (STAR21...) [HRP](#)  
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...) [DyLight®550](#), [DyLight®650](#), [DyLight®800](#)  
Goat Anti Rat IgG (STAR131...) [Alk. Phos.](#), [Biotin](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M402071:220718'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)