

Datasheet: 102005

BATCH NUMBER 157467

Description:	GOAT ANTI MOUSE IgM:HRP
Specificity:	IgM
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			1/4,000 - 1/8,000
Western Blotting			■	
Immunoblotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified Ig conjugated to Horseradish Peroxidase (HRP) - liquid
Antiserum Preparation	Purified Ig prepared by affinity chromatography on IgM covalently linked to agarose
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% Glycerol
Immunogen	Mouse IgM paraproteins
External Database Links	UniProt:

[P01872](#) [Related reagents](#)

[P01873](#) [Related reagents](#)

Entrez Gene:

[16019](#) Ighm [Related reagents](#)

[16019](#) Ighm [Related reagents](#)

RRID	AB_619903
Specificity	<p>Goat anti mouse IgM antibody recognises the heavy chain of mouse IgM as demonstrated by ELISA and flow cytometry. Minimal cross reactivity is observed with human immunoglobulins.</p> <p>Goat anti mouse IgM antibody has been cross absorbed against Mouse IgG₁, IgG_{2a}, IgG_{2b}, IgG₃ and IgA, pooled human sera and purified human paraproteins.</p>
References	<ol style="list-style-type: none">1. Kerr, K. <i>et al.</i> (2010) Inflammatory cytokine responses in a pregnant mouse model of <i>Chlamydomphila abortus</i> infection. Vet Microbiol. 144 (3-4): 392-8.2. Kamat, M.M. <i>et al.</i> (2016) Changes in Myeloid Lineage Cells in the Uterus and Peripheral Blood of Dairy Heifers During Early Pregnancy. Biol Reprod. 95 (3): 68.3. Monzo, H.J. <i>et al.</i> (2017) Insulin promotes cell migration by regulating PSA-NCAM. Exp Cell Res. 355 (1): 26-39.4. Desancé, M. <i>et al.</i> (2018) Chondrogenic Differentiation of Defined Equine Mesenchymal Stem Cells Derived from Umbilical Cord Blood for Use in Cartilage Repair Therapy. Int J Mol Sci. 19 (2)Feb 10 [Epub ahead of print].5. Penadés, M. <i>et al.</i> (2018) Long-term implications of feed energy source in different genetic types of reproductive rabbit females. II. Immunologic status. Animal. 12 (9): 1877-85.6. Penadés, M. <i>et al.</i> (2019) Early deviations in performance, metabolic and immunological indicators affect stayability in rabbit females. Animal. : 1-10.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	<p>Material Safety Datasheet documentation #10106 available at: https://www.bio-rad-antibodies.com/SDS/102005</p> <p>10106</p>
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)
[TMB CORE \(BUF056A\)](#)
[TMB CORE+ \(BUF062A\)](#)
[TMB SIGNAL+ \(BUF054A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350069:190307'

Printed on 25 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)