

Datasheet: 102005

**BATCH NUMBER 173009**

<b>Description:</b>	GOAT ANTI MOUSE IgM:HRP
<b>Specificity:</b>	IgM
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			1/4,000 - 1/8,000
Western Blotting			▪	
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified Ig conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Antiserum Preparation</b>	Purified Ig prepared by affinity chromatography on IgM covalently linked to agarose
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Mouse IgM paraproteins
<b>External Database Links</b>	<b>UniProt:</b>

[P01872](#) [Related reagents](#)

[P01873](#) [Related reagents](#)

**Entrez Gene:**

[16019](#) Ighm [Related reagents](#)

[16019](#) Ighm [Related reagents](#)

---

**RRID** AB\_619903

---

**Specificity** **Goat anti mouse IgM antibody** recognises the heavy chain of mouse IgM as demonstrated by ELISA and flow cytometry. Minimal cross reactivity is observed with human immunoglobulins.

Goat anti mouse IgM antibody has been cross absorbed against Mouse IgG<sub>1</sub>, IgG<sub>2a</sub>, IgG<sub>2b</sub>, IgG<sub>3</sub> and IgA, pooled human sera and purified human paraproteins.

---

**References**

1. Kerr, K. *et al.* (2010) Inflammatory cytokine responses in a pregnant mouse model of *Chlamydophila abortus* infection. [Vet Microbiol. 144 \(3-4\): 392-8.](#)
2. Kamat, M.M. *et al.* (2016) Changes in Myeloid Lineage Cells in the Uterus and Peripheral Blood of Dairy Heifers During Early Pregnancy. [Biol Reprod. 95 \(3\): 68.](#)
3. Monzo, H.J. *et al.* (2017) Insulin promotes cell migration by regulating PSA-NCAM. [Exp Cell Res. 355 \(1\): 26-39.](#)
4. Desancé, M. *et al.* (2018) Chondrogenic Differentiation of Defined Equine Mesenchymal Stem Cells Derived from Umbilical Cord Blood for Use in Cartilage Repair Therapy. [Int J Mol Sci. 19 \(2\): 537.](#)
5. Penadés, M. *et al.* (2018) Long-term implications of feed energy source in different genetic types of reproductive rabbit females. II. Immunologic status. [Animal. 12 \(9\): 1877-85.](#)
6. Penadés, M. *et al.* (2019) Early deviations in performance, metabolic and immunological indicators affect stayability in rabbit females. [Animal. : 1-10.](#)
7. Tchen, J. *et al.* (2024) PD-L1- and IL-4-expressing basophils promote pathogenic accumulation of T follicular helper cells in lupus. [Nat Commun. 15 \(1\): 3389.](#)

---

**Storage** Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

---

**Health And Safety Information** Material Safety Datasheet documentation #10106 available at: <https://www.bio-rad-antibodies.com/SDS/102005>

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M427877:240301'

**Printed on 29 Jan 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)