

Datasheet: 1021501

BATCH NUMBER 162199

| | |
|----------------------|-------------------------------|
| Description: | RAT ANTI MOUSE INTERLEUKIN-17 |
| Specificity: | IL-17 |
| Other names: | IL-17A |
| Format: | Purified |
| Product Type: | Monoclonal Antibody |
| Clone: | TC11-18H10 |
| Isotype: | IgG1 |
| Quantity: | 0.5 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | |
| ELISA | ▪ | | | |
| Western Blotting | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------------|-----------------------------|
| Target Species | Mouse |
| Product Form | Purified IgG - liquid |
| Buffer Solution | Borate buffered saline |
| Preservative Stabilisers | None present |
| Approx. Protein Concentrations | IgG concentration 0.5 mg/ml |
| RRID | AB_618680 |

Specificity **Rat anti Mouse interleukin-17 antibody, clone TC11-18H10** recognizes mouse interleukin 17, a ~32 kDa disulfide-linked homodimer consisting of two 136 amino acid chains, which belong to the IL-17 family (classified as IL-17A to IL-17F), and signal through the IL-17 receptor (IL-17R/CDw217). IL-17 is produced primarily by activated CD4⁺ T cells and induces a variety of target cells to secrete inflammatory cytokines, including IL-1beta, IL-6, TNF-alpha, IFN-gamma and granulocyte colony-stimulating factor, as well as synergizing with other inflammatory cytokines and agonists.

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10123 available at: <https://www.bio-rad-antibodies.com/SDS/1021501>
10123

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

| | |
|--|---|
| Rabbit Anti Rat IgG (STAR16...) | DyLight®800 |
| Rabbit Anti Rat IgG (STAR17...) | FITC |
| Goat Anti Rat IgG (STAR72...) | HRP |
| Goat Anti Rat IgG (STAR69...) | FITC |
| Goat Anti Rat IgG (STAR73...) | RPE |
| Rabbit Anti Rat IgG (STAR21...) | HRP |
| Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...) | DyLight®550 , DyLight®650 , DyLight®800 |
| Goat Anti Rat IgG (STAR131...) | Alk. Phos. , Biotin |

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M391999:211020'

Printed on 18 Jan 2024