

Datasheet: 104004

BATCH NUMBER 172450

Description:	GOAT ANTI MOUSE IgA:Alk. Phos.
Specificity:	IgA
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/2000 - 1/4000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Product Form

Purified IgG conjugated to Alkaline Phosphatase - liquid

Antiserum Preparation

Antiserum to mouse IgA was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

TRIS buffered saline, 1mM MgCl₂

Preservative Stabilisers

<0.1% Sodium Azide (NaN₃)
50% Glycerol

External Database Links

UniProt:

[P01878](#)

[Related reagents](#)

Entrez Gene:

[16061](#)

Igh-VJ558

[Related reagents](#)

RRID AB_619826

Specificity **Goat anti Mouse IgA antibody** recognizes murine IgA as demonstrated by ELISA, FLISA and/or Flow cytometry.
Goat anti Mouse IgA antibody has been cross adsorbed against mouse IgG1, IgG2a, IgG2b, IgG3 and IgM.

References

1. Lindner, C. *et al.* (2012) Age, microbiota, and T cells shape diverse individual IgA repertoires in the intestine. [J Exp Med. 209: 365-77.](#)
2. Tregoning, J.S. *et al.* (2013) A "prime-pull" vaccine strategy has a modest effect on local and systemic antibody responses to HIV gp140 in mice. [PLoS One. 8: e80559.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10322 available at: <https://www.bio-rad-antibodies.com/SDS/104004>

Regulatory For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M427890:240301'

Printed on 29 Jan 2026
