

Datasheet: 105005

BATCH NUMBER 163506

Description:	GOAT ANTI MOUSE KAPPA LIGHT CHAIN:HRP
Specificity:	KAPPA LIGHT CHAIN
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
ELISA	▪			1/4000 - 1/8000
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to mouse kappa light chain were made from pooled antisera from goats immunised with mouse paraproteins bearing kappa light chains. Antiserum was purified by affinity chromatography on pooled mouse Igs with kappa light chains. The antiserum was cross adsorbed on pooled mouse Igs with lambda light chains.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% glycerol.
Immunogen	Mouse paraproteins with kappa light chains.

**External Database
Links**

UniProt:

[P01837](#)

[Related reagents](#)

Entrez Gene:

[16071](#)

Igk-C

[Related reagents](#)

RRID

AB_609753

Specificity

Goat anti mouse kappa light chain polyclonal antibody is specific for mouse kappa light chains (free and bound). This antibody is cross-adsorbed against mouse lambda light chains.

References

1. Blakney, A.K. *et al.* (2021) Innate Inhibiting Proteins Enhance Expression and Immunogenicity of Self-Amplifying RNA. [Mol Ther. 29 \(3\): 1174-85.](#)

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

**Health And Safety
Information**

Material Safety Datasheet documentation #10106 available at:
<https://www.bio-rad-antibodies.com/SDS/105005>
10106

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402066:220718'

Printed on 01 Mar 2024