

Datasheet: 106001

BATCH NUMBER 158471

Description:	GOAT ANTI MOUSE LAMBDA LIGHT CHAIN
Specificity:	LAMBDA LIGHT CHAIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			
Immunofluorescence	▪			
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Mouse

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to mouse lambda light chains were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose.

Buffer Solution Borate buffered saline

Preservative Stabilisers None Present

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Mouse paraproteins with lambda 1 and lambda 2 light chains.

External Database Links

UniProt:

[P01843](#) [Related reagents](#)
[P01844](#) [Related reagents](#)

Entrez Gene:

[433053](#) LOC433053 [Related reagents](#)
[110786](#) Iglc2 [Related reagents](#)

RRID AB_609766

Specificity **Goat anti Mouse Lambda Light Chain antibody** recognizes mouse lambda 1 and lambda 2 light chains (free and bound) as demonstrated by ELISA. To minimise cross reactivity this product has been cross absorbed against mouse paraproteins with kappa light chains.

References

1. Rops, A.L. *et al.* (2010) The tetraspanin CD37 protects against glomerular IgA deposition and renal pathology. [Am J Pathol. 176 \(5\): 2188-97.](#)
2. Lewis, A. *et al.* (2017) Recovery from desensitization of IgE-dependent responses in human lung mast cells. [Clin Exp Allergy. Feb 25. \[Epub ahead of print\]](#)

Storage

Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10123 available at: <https://www.bio-rad-antibodies.com/SDS/106001>
10123

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350089:190307'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)