

Datasheet: 106001

BATCH NUMBER 161753

Description:	GOAT ANTI MOUSE LAMBDA LIGHT CHAIN
Specificity:	LAMBDA LIGHT CHAIN
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/4000 - 1/8000
Western Blotting	▪			
Immunofluorescence	▪			
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to mouse lambda light chains were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None Present

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Mouse paraproteins with lambda 1 and lambda 2 light chains.
External Database Links	<p>UniProt:</p> <p>P01843 Related reagents</p> <p>P01844 Related reagents</p> <p>Entrez Gene:</p> <p>433053 LOC433053 Related reagents</p> <p>110786 Iglc2 Related reagents</p>
RRID	AB_609766
Specificity	Goat anti Mouse Lambda Light Chain antibody recognizes mouse lambda 1 and lambda 2 light chains (free and bound) as demonstrated by ELISA. To minimise cross reactivity this product has been cross absorbed against mouse paraproteins with kappa light chains.
References	<ol style="list-style-type: none"> 1. Rops, A.L. <i>et al.</i> (2010) The tetraspanin CD37 protects against glomerular IgA deposition and renal pathology. Am J Pathol. 176 (5): 2188-97. 2. Lewis, A. <i>et al.</i> (2017) Recovery from desensitization of IgE-dependent responses in human lung mast cells. Clin Exp Allergy. Feb 25. [Epub ahead of print] 3. Blakney, A.K. <i>et al.</i> (2021) Innate Inhibiting Proteins Enhance Expression and Immunogenicity of Self-Amplifying RNA. Mol Ther. 29 (3): 1174-85.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10123 available at: https://www.bio-rad-antibodies.com/SDS/10600110123
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

'M392002:211020'

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)