

Datasheet: 106003

| | |
|----------------------|--|
| Description: | GOAT ANTI MOUSE LAMBDA LIGHT CHAIN:TRITC |
| Specificity: | LAMBDA LIGHT CHAIN |
| Format: | TRITC |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 0.5 mg |

Product Details

RRID AB_609770

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| Immunofluorescence | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Mouse

Product Form IgG conjugated to Rhodamine (TRITC) - liquid

| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
|-----------|-------------|---------------------|-------------------|
| | TRITC | 557 | 576 |

Antiserum Preparation Antisera to mouse lambda light chains were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 0.5 mg/ml

Immunogen Mouse paraproteins with lambda 1 and lambda 2 light chains.

**External Database
Links**

UniProt:

[P01843](#) [Related reagents](#)

[P01844](#) [Related reagents](#)

Entrez Gene:

[433053](#) LOC433053 [Related reagents](#)

[110786](#) Iglc2 [Related reagents](#)

Specificity

Goat anti Mouse Lambda Light Chain antibody recognizes mouse lambda 1 and lambda 2 light chains (free and bound) as demonstrated by ELISA. To minimise cross reactivity this product has been cross absorbed against mouse paraproteins with kappa light chains.

References

1. Rops, A.L. *et al.* (2010) The tetraspanin CD37 protects against glomerular IgA deposition and renal pathology. [Am J Pathol. 176 \(5\): 2188-97.](#)
2. Lewis, A. *et al.* (2017) Recovery from desensitization of IgE-dependent responses in human lung mast cells. [Clin Exp Allergy. Feb 25. \[Epub ahead of print\]](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

Please see label for expiry date

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M339161:181219'

Printed on 02 Jan 2019