

Datasheet: 107007

Description:	GOAT ANTI MOUSE IgG1:Texas Red®
Specificity:	IgG1
Format:	Texas Red®
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse		
Product Form	Purified IgG conjugated to Texas Red - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Texas Red	596	615

Antiserum Preparation Antisera to mouse IgG1 were raised by repeated immunisation of goats with highly purified antigen. Purified Ig was prepared from whole serum by affinity chromatography on Mouse IgG1 covalently linked to agarose.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	<0.1% sodium azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	Mouse IgG1 paraproteins.
External Database Links	<p>UniProt:</p> <p>P01868 Related reagents</p> <p>P01869 Related reagents</p> <p>Entrez Gene:</p> <p>16017 Ighg1 Related reagents</p> <p>16017 Ighg1 Related reagents</p>
Synonyms	Igh-4
RRID	AB_609718
Specificity	<p>Goat anti Mouse IgG1 antibody recognizes mouse immunoglobulin isotype G1.</p> <p>Goat anti Mouse IgG1 antibody has been adsorbed against mouse IgM, IgG2a, IgG2b, IgG3, IgA and human sera to minimise cross-reactivity.</p>
Flow Cytometry	Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl
References	<p>1. Jackson, S.J. <i>et al.</i> (2004) Cannabinoid-mediated neuroprotection following interferon-gamma treatment in a three-dimensional mouse brain aggregate cell culture. Eur J Neurosci. 20: 2267-75.</p> <p>2. Kushwaha, V. & Kaur, S. (2021) Cross-protective efficacy of immuno-stimulatory recombinant <i>Brugia malayi</i>. protein HSP60 against the <i>Leishmania donovani</i>. in BALB/c mice. Biologicals. 72: 18-26.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.</p>
Guarantee	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
Acknowledgements	Texas Red is a trademark of Molecular Probes, Inc.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/107007 10040
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M420310:230706'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)