

Datasheet: 109008

BATCH NUMBER 131014

Description:	GOAT ANTI MOUSE IgG2b:Biotin
Specificity:	IgG2b
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1/5000 - 1/10000
Immunoprecipitation			■	
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to mouse IgG2b were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Mouse IgG2b paraproteins
External Database Links	<p>UniProt: P01867 Related reagents</p> <p>Entrez Gene: 16016 Ighg2b Related reagents</p>
RRID	AB_609729
Specificity	<p>Goat anti mouse IgG2b polyclonal antibody is specific for the heavy chain of mouse IgG2b.</p> <p>This antibody has been cross-adsorbed against mouse IgM, IgG1, IgG2a, IgG3, IgA and human serum to minimise cross-reactivity.</p>
Flow Cytometry	Use 50ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/109008 10040</p>
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M362900:200528'

Printed on 01 Mar 2024