

Datasheet: 12004165

Description:	ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY
Specificity:	TUBULIN
Format:	Rhodamine
Product Type:	Monoclonal Antibody
Clone:	AbD22584
Isotype:	HuCAL Fab monovalent
Quantity:	200 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1,000 - 1/10,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

A monovalent human recombinant Fab selected from the HuCAL® phage display library, expressed in *E. coli*. and purified using NiNTA affinity chromatography. This antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain and conjugated to rhodamine - lyophilized

Reconstitution

Resuspend contents of the tube in the indicated volume of distilled or deionized water and leave on ice for at least 30 min prior to use. Brief centrifugation (pulse spin for 2-3 sec at maximum speed on a tabletop microcentrifuge) may be used to collect the contents at the bottom of the tube.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
-----------	-------------	---------------------	-------------------

Buffer Solution Phosphate buffered saline

Preservative >80% Sucrose
Stabilisers <0.04% 2-Methyl-2H-Isothiazol-3-One
<10% Bovine Serum Albumin

External Database Links

UniProt:
[P07437](#) [Related reagents](#)

Entrez Gene:
[203068](#) TUBB [Related reagents](#)

Synonyms TUBB5

Product Information hFAB Rhodamine Housekeeping Protein Fluorescent Primary Antibodies allow one-step detection of common housekeeping proteins like actin, tubulin, and GAPDH in your samples. These anti-[actin](#), anti-tubulin, and anti-[GAPDH](#) antibodies are directly labeled with a rhodamine derivative (Exmax/Emmax = 530 nm/580 nm). The hFAB Rhodamine Antibodies were generated with [Human Combinatorial Antibody Library \(HuCAL®\) Technology](#), providing exceptional specificity.

These housekeeping protein antibodies can be used for protein normalization when the protein loading amounts are in the linear dynamic range. hFAB Rhodamine Antibodies can be imaged with the ChemiDoc™ MP Imaging System, and can be used in multiplex experiments with StarBright™ Blue 700 Secondary Antibodies as well as other traditional fluorescent antibodies, permitting detection of multiple target proteins and the housekeeping protein in a single blot.

Features and Benefits

Simple — No requirement for a secondary antibody; one-step detection of housekeeping proteins using highly cross-adsorbed antibodies for high specificity and sensitivity

Versatile — Compatible with use of a wide range of other fluorophores including StarBright Blue 700 and DyLight 800 Fluorophores for multiplex western blotting

Unique — HuCAL Technology eliminates cross-reactivity with primary/secondary antibodies from other species

Diverse — Antibodies have been tested to detect human or mouse housekeeping proteins

Instructions For Use Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/12004163.pdf

Storage 1 year at -20°C lyophilized; 6 months at +4°C after resuspension.
This product is photosensitive and should be protected from light
DO NOT FREEZE the solubilized material

Guarantee Guaranteed until date of expiry. Please see product label.

Acknowledgements Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
His-tag is a registered trademark of EMD Biosciences.

Health And Safety Information Material Safety Datasheet documentation #12004165 available at:
12004165: <https://www.bio-rad-antibodies.com/uploads/MSDS/12004165.pdf>
Material Safety Datasheet documentation #12004165 available at
<https://www.bio-rad-antibodies.com/uploads/MSDS/12004165.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[SHEEP ANTI RABBIT IgG:DyLight®800 \(STAR36D800GA\)](#)

[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004168\)](#)

[GOAT ANTI MOUSE IgG \(H/L\):DyLight®800 \(MULTI SPECIES ADSORBED\) \(STAR117D800GA\)](#)

[GOAT ANTI-MOUSE IgG StarBright™ Blue 700 \(12004158\)](#)

[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004161\)](#)

[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004167\)](#)

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004163\)](#)

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004164\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

'M374498:201109'

Printed on 02 Dec 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)