

Datasheet: 12004167

Description:	ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY
Specificity:	GAPDH
Format:	Rhodamine
Product Type:	Monoclonal Antibody
Clone:	AbD22549
Isotype:	HuCAL Fab monovalent
Quantity:	200 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1,000 - 1/10,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

A monovalent human recombinant Fab selected from the HuCAL® phage display library, expressed in *E. coli*. and purified using NiNTA affinity chromatography. This antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain and conjugated to rhodamine - lyophilized

Reconstitution

Resuspend contents of the tube in the indicated volume of distilled or deionized water and leave on ice for at least 30 min prior to use. Brief centrifugation (pulse spin for 2-3 sec at maximum speed on a tabletop microcentrifuge) may be used to collect the contents at the bottom of the tube.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Rhodamine	530	580
Buffer Solution	Phosphate buffered saline		
Preservative	>80% Sucrose		
Stabilisers	<0.04% 2-Methyl-2H-Isothiazol-3-One		
	<10% Bovine Serum Albumin		
External Database Links	UniProt: P04406 Related reagents		
	Entrez Gene: 2597 GAPDH Related reagents		
Synonyms	GAPD		
Product Information	<p>hFAB Rhodamine Housekeeping Protein Fluorescent Primary Antibodies allow one-step detection of common housekeeping proteins like actin, tubulin, and GAPDH in your samples. These anti-actin, anti-tubulin, and anti-GAPDH antibodies are directly labeled with a rhodamine derivative (Exmax/Emmax = 530 nm/580 nm). The hFAB Rhodamine Antibodies were generated with Human Combinatorial Antibody Library (HuCAL®) Technology, providing exceptional specificity.</p> <p>These housekeeping protein antibodies can be used for protein normalization when the protein loading amounts are in the linear dynamic range. hFAB Rhodamine Antibodies can be imaged with the ChemiDoc™ MP Imaging System, and can be used in multiplex experiments with StarBright™ Blue 700 Secondary Antibodies as well as other traditional fluorescent antibodies, permitting detection of multiple target proteins and the housekeeping protein in a single blot.</p> <p>Features and Benefits</p> <p>Simple — No requirement for a secondary antibody; one-step detection of housekeeping proteins using highly cross-adsorbed antibodies for high specificity and sensitivity</p> <p>Versatile — Compatible with use of a wide range of other fluorophores including StarBright Blue 700 and DyLight 800 Fluorophores for multiplex western blotting</p> <p>Unique — HuCAL Technology eliminates cross-reactivity with primary/secondary antibodies from other species</p> <p>Diverse — Antibodies have been tested to detect human or mouse housekeeping proteins</p>		
Instructions For Use	Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/12004163.pdf		
References	1. Meissner, M.E. <i>et al.</i> (2022) Differential Activity of APOBEC3F, APOBEC3G, and APOBEC3H in the Restriction of HIV-2. J Mol Biol. 434 (2): 167355.		

Storage	1 year at -20°C lyophilized; 6 months at +4°C after resuspension. This product is photosensitive and should be protected from light DO NOT FREEZE the solubilized material
Guarantee	Guaranteed until date of expiry. Please see product label.
Acknowledgements	This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #12004167 available at: 12004167: https://www.bio-rad-antibodies.com/uploads/MSDS/12004167.pdf Material Safety Datasheet documentation #12004167 available at https://www.bio-rad-antibodies.com/uploads/MSDS/12004167.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004163\)](#)

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004164\)](#)

[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004165\)](#)

[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004166\)](#)

[GOAT ANTI-MOUSE IgG StarBright™ Blue 700 \(12004158\)](#)

[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004161\)](#)

[GOAT ANTI MOUSE IgG \(H/L\):DyLight®800 \(MULTI SPECIES ADSORBED\) \(STAR117D800GA\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M391554:211015'

Printed on 21 Mar 2022