

Datasheet: 12024057

Description:	TrailBlazer SBY665 Label Kit
Name:	TrailBlazer StarBright Yellow 665 Label Kit
Format:	StarBright Yellow 665
Product Type:	Conjugation Kit
Quantity:	1 X 100 µg Antibody Conjugation

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Conjugation - Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	StarBright Yellow 665	554	670

Product Information

TrailBlazer™ StarBright™ Dye Label is the second part of a two-kit antibody labeling system. Together with the [TrailBlazer Tag kit](#) it will allow you to conjugate an antibody to superior StarBright Dyes using the innovative [SpyTag and SpyCatcher technology](#). They provide the benefits of using StarBright Dyes in flow cytometry and western blots for your antibody of choice.

First, the SpyTag is attached to an antibody using the TrailBlazer Tag Kit. This antibody, now with a SpyTag, is then mixed with a StarBright Dye labeled SpyCatcher using the Label Kit. A spontaneous reaction occurs between the SpyTag and SpyCatcher resulting in an irreversible covalent isopeptide bond. The result is a highly stable StarBright Dye labeled antibody, which should be stored at +4°C in the dark, and when stored correctly will be stable for up to 12 months.

Batch specific [Certificate of Analysis](#) for TrailBlazer StarBright Dye Label Kits are available to be downloaded directly from the website.

In flow cytometry we suggest using 5µl of the undiluted labeled antibody to stain 10⁶ cells in 100µl. However, we highly recommend that you [titrate the antibody](#) to determine the optimal amount of antibody for your specific experiment.

The TrailBlazer StarBright Dye Label Kits are available in different pack sizes and for a range of StarBright Dyes. Discover the [full range](#) and all the [benefits of these TrailBlazer Tag and StarBright Dye Label Kits](#).

Reagents In The Kit	1 X 100 µg Antibody Conjugation (3 X 100 µg Antibody Conjugation Kit available) 1 vial Storage Buffer 1 vial StarBright Dye-Catcher SBY665 1 vial Quenching Reagent 1 Spin Column 1 Collection Tube
----------------------------	--

Flow Cytometry	Use 5ul of the suggested working dilution to label 10 ⁶ cells in 100ul. Best practices suggest a 5 min centrifugation at 6,000g prior to sample application.
-----------------------	---

Instructions For Use

Add the SpyTagged antibody product from TrailBlazer Tag Kit to the StarBright Dye-Catcher vial. Mix and incubate for 30 min at room temperature or overnight at +4°C. Add the sample to the Spin Column. Vortex briefly and mix on an end-over-end rotary mixer at room temperature for 60 min.

Break open the bottom of the Spin Column and insert the column into the 2 ml amber collection tube. Slightly loosen the cap of the Spin Column and centrifuge at 1500 x g for 2 min. Retain the eluate, containing the label antibody.

Prepare Quenching Reagent by adding 5 µl of ultrapure water.

Add 4 µl of prepared Quenching Reagent and 45 µl of Storage Buffer to the labeled antibody. Vortex briefly to mix.

Antibodies are now labeled and ready to use. Store at +4°C in the dark for up to 1 year. For full instructions for use please go to: bio-rad-antibodies.com/label-SB-kit-IFU.

Storage	Store at +4°C. DO NOT FREEZE.
----------------	-------------------------------

The StarBright Dye labeled antibody should be stored undiluted at +4°C in the dark for up to 1 year.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	This product is covered by U.S. Patent No. 10,150,841 and related U.S. and foreign counterparts This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #12024057 available at: https://www.bio-rad-antibodies.com/SDS/12024057
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[TrailBlazer Tag Kit 0.1 mg \(12020038\)](#)

[TrailBlazer SBB700 Label Kit \(12020040\)](#)

[TrailBlazer SBV515 Label Kit \(12020041\)](#)

[TrailBlazer SBV440 Label Kit \(12020042\)](#)

[TrailBlazer SBV610 Label Kit \(12020043\)](#)
[TrailBlazer SBV670 Label Kit \(12020044\)](#)
[TrailBlazer SBV570 Label Kit \(12020045\)](#)
[TrailBlazer SBY605 Label Kit \(12023837\)](#)
[TrailBlazer SBB675 Label Kit \(12023949\)](#)
[TrailBlazer SBB580 Label Kit \(12023991\)](#)
[TrailBlazer SBB765 Label Kit \(12023960\)](#)
[TrailBlazer SBY800 Label Kit \(12024044\)](#)
[TrailBlazer SBY720 Label Kit \(12023838\)](#)
[TrailBlazer SBY575 Label Kit \(12023974\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M444665:250815'

Printed on 16 Dec 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)