

Datasheet: 17010822

Description:	TrailBlazer SBR815 Label Kit 3x 0.1 mg
Name:	TrailBlazer StarBright Red 815 Label Kit
Format:	StarBright Red 815
Product Type:	Conjugation Kit
Quantity:	3 X 100 µg Antibody Conjugations

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Conjugation - Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	StarBright Red 815	654	811

Product Information

TrailBlazer™ StarBright™ Dye Label is the second part of a two-kit antibody labeling system. Together with the [TrailBlazer Tag kit](#) it will allow you to conjugate an antibody to superior StarBright Dyes using the innovative [SpyTag and SpyCatcher technology](#). They provide the benefits of using StarBright Dyes in flow cytometry and western blots for your antibody of choice.

First, the SpyTag is attached to an antibody using the TrailBlazer Tag Kit. This antibody, now with a SpyTag, is then mixed with a StarBright Dye labeled SpyCatcher using the Label Kit. A spontaneous reaction occurs between the SpyTag and SpyCatcher resulting in an irreversible covalent isopeptide bond. The result is a highly stable StarBright Dye labeled antibody, which should be stored at +4°C in the dark, and when stored correctly will be stable for up to 12 months.

Batch specific [Certificate of Analysis](#) for TrailBlazer StarBright Dye Label Kits are available to be downloaded directly from the website.

In flow cytometry we suggest using 5µl of the undiluted labeled antibody to stain 10⁶ cells in 100µl. However, we highly recommend that you [titrate the antibody](#) to determine the optimal amount of antibody for your specific experiment.

The TrailBlazer StarBright Dye Label Kits are available in different pack sizes and for a range of StarBright Dyes. Discover the [full range](#) and all the [benefits of these TrailBlazer Tag and StarBright Dye Label Kits](#).

Reagents In The Kit 3 X 100 µg Antibody Conjugations ([1 X 100 µg Antibody Conjugation Kit available](#))
3 vials Storage Buffer
3 vials StarBright Dye-Catcher SBR815
3 vials Quenching Reagent
3 Spin Columns
3 Collection Tubes

Flow Cytometry Use 5ul of the suggested working dilution to label 10⁶ cells in 100ul. Best practices suggest a 5 min centrifugation at 6,000g prior to sample application.

Instructions For Use

Add the SpyTagged antibody product from TrailBlazer Tag Kit to the StarBright Dye-Catcher vial. Mix and incubate for 30 min at room temperature or overnight at +4°C.

Add the sample to the Spin Column. Vortex briefly and mix on an end-over-end rotary mixer at room temperature for 60 min.

Break open the bottom of the Spin Column and insert the column into the 2 ml amber collection tube. Slightly loosen the cap of the Spin Column and centrifuge at 1500 x g for 2 min. Retain the eluate, containing the label antibody.

Prepare Quenching Reagent by adding 5 µl of ultrapure water.

Add 4 µl of prepared Quenching Reagent and 45 µl of Storage Buffer to the labeled antibody. Vortex briefly to mix.

Antibodies are now labeled and ready to use. Store at +4°C in the dark for up to 1 year. For full instructions for use please go to: bio-rad-antibodies.com/label-SB-kit-IFU.

Storage Store at +4°C. DO NOT FREEZE.

The StarBright Dye labeled antibody should be stored undiluted at +4°C in the dark for up to 1 year.

Guarantee 12 months from date of despatch

Acknowledgements This product is covered by U.S. Patent No. 10,150,841 and related U.S. and foreign counterparts This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details

Health And Safety Information Material Safety Datasheet documentation #17010822 available at: <https://www.bio-rad-antibodies.com/SDS/17010822>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[TrailBlazer Tag Kit 0.1 mg \(12020038\)](#)

[TrailBlazer Tag Kit 3 x 0.1 mg \(12020901\)](#)

[TrailBlazer SBR815 Label Kit \(12024063\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M446194:251022'

Printed on 06 Jan 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)