

Datasheet: 1707-0189G

BATCH NUMBER 1608

| | |
|----------------------|------------------------------------|
| Description: | GOAT ANTI HUMAN C-REACTIVE PROTEIN |
| Specificity: | C-REACTIVE PROTEIN |
| Other names: | CRP |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-----------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | 1/1000 |
| Immunodiffusion | ▪ | | | 1/30 - 1/200 |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Human

Product Form Purified IgG - liquid

Antiserum Preparation Antisera to human C-reactive protein were raised by repeated immunisations of goats with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 5.0mg/ml

Immunogen Native human C-reactive protein.

**External Database
Links**

UniProt:

[P02741](#) [Related reagents](#)

Entrez Gene:

[1401](#) CRP [Related reagents](#)

Synonyms

PTX1

RRID

AB_2057775

Specificity

Goat anti Human C-reactive protein antibody recognizes human C-reactive protein (CRP), a calcium-binding acute phase reactant belonging to the pentraxin family which is secreted into the plasma by the liver, in response to cytokine stimulation.

In addition to binding to phosphocholine residues, CRP can bind to other autologous and extrinsic ligands on targets, including damaged cell membranes, ribonucleoprotein particles, apoptotic cells and constituents of micro-organisms, and may activate the classical complement pathway through C1q recognition, making CRP a factor in a variety of host defense mechanisms. Levels of plasma CRP are significantly increased during the acute phase response, tissue injury and infection and can be used diagnostically as a marker of inflammation, disease progress and risk assessment of cardiovascular disease.

References

1. Günter, G. *et al.* (2019) Smartphone biosensor for Salmonella and Amitriptyline [Proc. SPIE 11028, Optical Sensors 2019, 110280R](#)

Further Reading

1. Pepys, M.B. & Hirschfield, G.M. (2003) C-reactive protein: a critical update. [J Clin Invest. 111 \(12\): 1805-12.](#)
2. Clearfield, M.B. (2005) C-reactive protein: a new risk assessment tool for cardiovascular disease. [J Am Osteopath Assoc. 105 \(9\): 409-16.](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M362921:200528'

Printed on 08 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)