

Datasheet: 1707-2029

BATCH NUMBER 162651

Description:	NATIVE HUMAN C-REACTIVE PROTEIN
Name:	C-REACTIVE PROTEIN
Other names:	CRP
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified C-reactive Protein - liquid
Preparation	Purified protein prepared by affinity chromatography on CH-Sepharose and hydroxylapatite
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Total protein concentration 1.0 mg/ml

External Database Links

UniProt:
[P02741](https://www.uniprot.org/entry/P02741) [Related reagents](#)

Entrez Gene:

[1401](#) CRP [Related reagents](#)

Synonyms	PTX1
-----------------	------

Product Information	Human C-reactive protein is a sterile filtered preparation derived from pleural ascites. C-reactive protein is an acute phase protein produced by the liver, a marker of inflammation and useful for monitoring and predicting coronary artery disease.
----------------------------	--

Protein Molecular Weight	May exist as a 140 kDa pentamer under native conditions or as a 26 kDa monomer under reduced conditions
---------------------------------	---

Purity	>98% by SDS PAGE
---------------	------------------

References	<ol style="list-style-type: none">1. Tomita, S. <i>et al.</i> (2016) Artificial Modification of an Enzyme for Construction of Cross-Reactive Polyion Complexes To Fingerprint Signatures of Proteins and Mammalian Cells. Anal Chem. 88 (18): 9079-86.2. Zheng, Y. & Demarco M. (2017) Manipulating trypsin digestion conditions to accelerate proteolysis and simplify digestion workflows in development of protein mass spectrometric assays for the clinical laboratory Clinical Mass Spectrometry. [Epub ahead of print].
-------------------	---

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
----------------	---

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10327 available at: https://www.bio-rad-antibodies.com/SDS/1707-2029 10327 Donor material tested and found negative for HBsAg, HCV, and HIV1 and 2 antibodies. As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350147:190307'

Printed on 20 Mar 2024