

Datasheet: 1707-2029

**BATCH NUMBER 167519**

<b>Description:</b>	NATIVE HUMAN C-REACTIVE PROTEIN
<b>Name:</b>	C-REACTIVE PROTEIN
<b>Other names:</b>	CRP
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified C-reactive Protein - liquid
<b>Preparation</b>	Purified protein prepared by affinity chromatography on CH-Sepharose and hydroxylapatite
<b>Buffer Solution</b>	TRIS buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Total protein concentration 1.0 mg/ml

### External Database Links

**UniProt:**  
[P02741](https://www.uniprot.org/entry/P02741)    [Related reagents](#)

**Entrez Gene:**

[1401](#) CRP [Related reagents](#)

---

<b>Synonyms</b>	PTX1
-----------------	------

---

<b>Product Information</b>	<b>Human C-reactive protein</b> is a sterile filtered preparation derived from pleural ascites. C-reactive protein is an acute phase protein produced by the liver, a marker of inflammation and useful for monitoring and predicting coronary artery disease.
----------------------------	--

---

<b>Protein Molecular Weight</b>	May exist as a 140 kDa pentamer under native conditions or as a 26 kDa monomer under reduced conditions
---------------------------------	---

---

<b>Purity</b>	>98% by SDS PAGE
---------------	------------------

---

<b>References</b>	<ol style="list-style-type: none"><li>Tomita, S. <i>et al.</i> (2016) Artificial Modification of an Enzyme for Construction of Cross-Reactive Polyion Complexes To Fingerprint Signatures of Proteins and Mammalian Cells. <a href="#">Anal Chem. 88 (18): 9079-86.</a></li><li>Zheng, Y. &amp; Demarco M. (2017) Manipulating trypsin digestion conditions to accelerate proteolysis and simplify digestion workflows in development of protein mass spectrometric assays for the clinical laboratory <a href="#">Clinical Mass Spectrometry. [Epub ahead of print].</a></li></ol>
-------------------	---

---

<b>Storage</b>	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
----------------	---

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10327 available at: <a href="https://www.bio-rad-antibodies.com/SDS/1707-2029">https://www.bio-rad-antibodies.com/SDS/1707-2029</a> 10327  Donor material tested and found negative for HBsAg, HCV, and HIV1 and 2 antibodies.  As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M426255:231207'

Printed on 20 Mar 2024