

Datasheet: 201001

BATCH NUMBER 166074

Description:	GOAT ANTI HUMAN Ig
Specificity:	Ig
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Western Blotting			■	
Immunoblotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to Human Ig was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	1.0 mg/ml
Immunogen	Normal human IgM and IgG, and human IgM, IgG and IgA paraproteins.

Specificity	<p>Goat anti human Ig polyclonal antibody reacts with the heavy and light chains of human IgM, IgG and IgA, and with the light chains of human IgD as demonstrated by ELISA and flow cytometry.</p> <p>This antibody may also react with immunoglobulins from other species</p>
References	<p>1. Lavender, K.J. <i>et al.</i> (2013) BLT-humanized C57BL/6 Rag2-/-yc-/-CD47-/- mice are resistant to GVHD and develop B- and T-cell immunity to HIV infection. Blood. 122: 4013-20.</p>
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	<p>Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.</p>
Health And Safety Information	<p>Material Safety Datasheet documentation #10123 available at: https://www.bio-rad-antibodies.com/SDS/20100110123</p>
Regulatory	<p>For research purposes only</p>

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America	<p>Tel: +1 800 265 7376</p> <p>Fax: +1 919 878 3751</p> <p>Email: antibody_sales_us@bio-rad.com</p>	Worldwide	<p>Tel: +44 (0)1865 852 700</p> <p>Fax: +44 (0)1865 852 739</p> <p>Email: antibody_sales_uk@bio-rad.com</p>	Europe	<p>Tel: +49 (0) 89 8090 95 21</p> <p>Fax: +49 (0) 89 8090 95 50</p> <p>Email: antibody_sales_de@bio-rad.com</p>
----------------------------------	--	------------------	--	---------------	--

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M401882:220718'

Printed on 01 Mar 2024