

Datasheet: 202004

Description:	GOAT ANTI HUMAN IgM:Alk. Phos.
Specificity:	IgM
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			1/2000 - 1/4000
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Alkaline Phosphatase - liquid
Antiserum Preparation	Antiserum to human IgM was raised by repeated immunisation of goats with human IgM paraproteins. The preparation was purified by affinity chromatography on human IgM covalently linked to agarose.
Buffer Solution	TRIS buffered saline, 1mM MgCl ₂
Preservative	0.1% Sodium Azide
Stabilisers	50% Glycerol
Immunogen	Human IgM
External Database Links	<p>UniProt: P01871 Related reagents</p> <p>Entrez Gene: 3507 IGHM Related reagents</p>
RRID	AB_619808

Specificity	Goat anti human IgM polyclonal antibody reacts with the heavy chain of human IgM. This antibody has been cross-absorbed to reduce cross-reactivity with human IgG and IgA, but may react with IgM from other species.
--------------------	--

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10322 available at: 10322: https://www.bio-rad-antibodies.com/uploads/MSDS/10322.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M360335:191106'

Printed on 01 May 2020