

Datasheet: 202009

**BATCH NUMBER 157074**

<b>Description:</b>	GOAT ANTI HUMAN IgM:RPE
<b>Specificity:</b>	IgM
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578

**Antiserum Preparation** Antiserum to human IgM was raised by repeated immunisation of goats with human IgM paraproteins. The preparation was purified by affinity chromatography on human IgM covalently linked to Sepharose-4B.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.1% Sodium Azide
<b>Stabilisers</b>	Stabilising agent
<b>Approx. Protein Concentrations</b>	Ig concentration 0.5 mg/ml

Immunogen Human IgM

---

**External Database**

**Links**

**UniProt:**

[P01871](#) [Related reagents](#)

**Entrez Gene:**

[3507](#) IGHM [Related reagents](#)

---

RRID AB\_616853

---

**Specificity**

**Goat anti human IgM polyclonal antibody** reacts with the heavy chain of human IgM. This antibody has been cross-absorbed to reduce cross-reactivity with human IgG and IgA, but may react with IgM from other species.

---

**Flow Cytometry**

Use 10ul of the suggested working dilution to label  $1 \times 10^6$  cells in 100ul.

---

**Storage**

Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10331 available at:  
<https://www.bio-rad-antibodies.com/SDS/20200910331>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350154:190307'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)