

Datasheet: 204005

**BATCH NUMBER 155479**

<b>Description:</b>	GOAT ANTI HUMAN IgG:HRP
<b>Specificity:</b>	IgG
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
ELISA	▪			1/4,000 - 1/8,000
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Human

### Product Form

Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

### Antiserum Preparation

Antisera to human IgG were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

50% Glycerol

### Approx. Protein Concentrations

IgG concentration 0.6 mg/ml

### Immunogen

Normal human IgG and human IgG paraproteins

**External Database  
Links**

**UniProt:**

[P01857](#)    [Related reagents](#)  
[P01859](#)    [Related reagents](#)  
[P01861](#)    [Related reagents](#)  
[P01860](#)    [Related reagents](#)

**Entrez Gene:**

[3500](#)    IGHG1    [Related reagents](#)  
[3501](#)    IGHG2    [Related reagents](#)  
[3502](#)    IGHG3    [Related reagents](#)  
[3503](#)    IGHG4    [Related reagents](#)

---

**RRID**                      AB\_619881

---

**Specificity**              **Goat anti Human IgG antibody** reacts with the heavy chain of human IgG and has been cross-absorbed with human IgM and IgA, but may recognise IgG from other species.

---

**References**              1. Lund K *et al.* (2014) The novel EpCAM-targeting monoclonal antibody 3-17I linked to saporin is highly cytotoxic after photochemical internalization in breast, pancreas and colon cancer cell lines. [MAbs. 6 \(4\): 1038-50.](#)  
2. Gilbert, A. (2013) Humoral immune response to melanoma: discovery and evaluation of anti-melanoma antibodies [PhD. Thesis KCL UK.](#)

---

**Storage**                      Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**                      Guaranteed until date of expiry. Please see product label.

---

**Health And Safety  
Information**                  Material Safety Datasheet documentation #10106 available at:  
<https://www.bio-rad-antibodies.com/SDS/204005>  
10106

---

**Regulatory**                      For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)  
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)  
[TMB CORE \(BUF056A\)](#)  
[TMB CORE+ \(BUF062A\)](#)  
[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**Worldwide** Tel: +44 (0)1865 852 700

**Europe** Tel: +49 (0) 89 8090 95 21

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M360671:200109'

**Printed on 17 Apr 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)