

Datasheet: 204009

**BATCH NUMBER 159724**

<b>Description:</b>	GOAT ANTI HUMAN IgG:RPE
<b>Specificity:</b>	IgG
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578
<b>Preparation</b>	Purified IgG prepared by affinity chromatography		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> ) Stabilizing agent (sucrose)		
<b>Approx. Protein Concentrations</b>	0.5 mg/ml		
<b>Immunogen</b>	Normal human IgG and human IgG paraproteins		

**External Database  
Links**

**UniProt:**

[P01857](#)    [Related reagents](#)  
[P01859](#)    [Related reagents](#)  
[P01861](#)    [Related reagents](#)  
[P01860](#)    [Related reagents](#)

**Entrez Gene:**

[3500](#) IGHG1    [Related reagents](#)  
[3501](#) IGHG2    [Related reagents](#)  
[3502](#) IGHG3    [Related reagents](#)  
[3503](#) IGHG4    [Related reagents](#)

---

**RRID**                      AB\_616852

---

**Specificity**              **Goat anti Human IgG antibody** reacts with the heavy chain of human IgG and has been cross-absorbed with human IgM and IgA, but may recognise IgG from other species.

---

**References**              1. Lund K *et al.* (2014) The novel EpCAM-targeting monoclonal antibody 3-17I linked to saporin is highly cytotoxic after photochemical internalization in breast, pancreas and colon cancer cell lines. [MAbs. 6 \(4\): 1038-50.](#)  
2. Gilbert, A. (2013) Humoral immune response to melanoma: discovery and evaluation of anti-melanoma antibodies [PhD. Thesis KCL UK.](#)

---

**Storage**                      Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**                      Guaranteed until date of expiry. Please see product label.

---

**Health And Safety  
Information**              Material Safety Datasheet documentation #10045 available at:  
<https://www.bio-rad-antibodies.com/SDS/204009>  
10045

---

**Regulatory**                      For research purposes only

---

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)  
[HUMAN SEROBLOCK \(BUF070B\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
--------------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M387559:210629'

