

Datasheet: 204009

**BATCH NUMBER 173676**

<b>Description:</b>	GOAT ANTI HUMAN IgG:RPE
<b>Specificity:</b>	IgG
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578
<b>Preparation</b>	Purified IgG prepared by affinity chromatography		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> ) Stabilizing agent (sucrose)		
<b>Approx. Protein Concentrations</b>	0.5 mg/ml		
<b>Immunogen</b>	Normal human IgG and human IgG paraproteins		

**External Database  
Links**

**UniProt:**

[P01857](#)    [Related reagents](#)  
[P01859](#)    [Related reagents](#)  
[P01861](#)    [Related reagents](#)  
[P01860](#)    [Related reagents](#)

**Entrez Gene:**

[3500](#)    IGHG1    [Related reagents](#)  
[3501](#)    IGHG2    [Related reagents](#)  
[3502](#)    IGHG3    [Related reagents](#)  
[3503](#)    IGHG4    [Related reagents](#)

---

**RRID**                      AB\_616852

---

**Specificity**              **Goat anti Human IgG antibody** reacts with the heavy chain of human IgG and has been cross-absorbed with human IgM and IgA, but may recognise IgG from other species.

---

**References**              1. Lund K *et al.* (2014) The novel EpCAM-targeting monoclonal antibody 3-17I linked to saporin is highly cytotoxic after photochemical internalization in breast, pancreas and colon cancer cell lines. [MAbs. 6 \(4\): 1038-50.](#)  
2. Gilbert, A. (2013) Humoral immune response to melanoma: discovery and evaluation of anti-melanoma antibodies [PhD. Thesis KCL UK.](#)

---

**Storage**                      This product is shipped at ambient temperature.  
Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**                      Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

---

**Health And Safety  
Information**                      Material Safety Datasheet documentation #10045 available at:  
<https://www.bio-rad-antibodies.com/SDS/204009>

---

**Regulatory**                      For research purposes only

---

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)  
[HUMAN SEROBLOCK \(BUF070B\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M438842:250523'

