

Datasheet: 205004

BATCH NUMBER 165121

Description:	GOAT ANTI HUMAN IgA:Alk. Phos.
Specificity:	IgA
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/2000 - 1/4000
Immunoprecipitation			▪	
Western Blotting			▪	
Immunoblotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Alkaline Phosphatase - liquid
Preparation	Purified IgG prepared by affinity chromatography
Antiserum Preparation	Antisera to human IgA were raised by repeated immunisations of goats with highly purified antigen.
Buffer Solution	TRIS buffered saline, 1mM MgCl ₂
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃) 50% Glycerol

Immunogen	Human IgA.
External Database Links	<p>UniProt:</p> <p>P01876 Related reagents</p> <p>P01877 Related reagents</p> <p>Entrez Gene:</p> <p>3493 IGHA1 Related reagents</p> <p>3494 IGHA2 Related reagents</p>
RRID	AB_609670
Specificity	<p>Goat anti Human IgA antibody recognizes the heavy chain of human IgA as demonstrated by ELISA and Flow cytometry. IgA is the major category of antibody present in mucosal secretions and acts as a first line of defense against inhaled or ingested pathogens at the mucosal interface. IgA is also present in serum where it acts to eliminate pathogens via interaction with FcαR, also known as CD89. (Woof et al. 2004). Biotinylated Goat anti Human IgA antibody (205008) has been cross-adsorbed against human IgD, IgG and IgM to ensure class specificity.</p>
References	<ol style="list-style-type: none"> 1. Wright, A.K. et al. (2012) Human Nasal Challenge with <i>Streptococcus pneumoniae</i> Is Immunising in the Absence of Carriage. PLoS Pathog. 8: e1002622. 2. Pennington, S.H. et al. (2016) Polysaccharide-specific Memory B-cells Predict Protection Against Experimental Human Pneumococcal Carriage. Am J Respir Crit Care Med. Jul 12. [Epub ahead of print] 3. Turner, P.J. et al. (2020) Differences in nasal immunoglobulin A responses to influenza vaccine strains after live attenuated influenza vaccine (LAIV) immunization in children. Clin Exp Immunol. 199 (2): 109-118. 4. Cole, M.E. et al. (2021) Pre-existing influenza-specific nasal IgA or nasal viral infection does not affect live attenuated influenza vaccine immunogenicity in children. Clin Exp Immunol. 204 (1): 125-133.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10322 available at: https://www.bio-rad-antibodies.com/SDS/205004 10322
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M402043:220718'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)