

Datasheet: 205009

**BATCH NUMBER 158373**

<b>Description:</b>	GOAT ANTI HUMAN IgA:RPE
<b>Specificity:</b>	IgA
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on pooled human IgA		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative</b>	0.1% Sodium azide		
<b>Stabilisers</b>	Stabilising agent		
<b>Approx. Protein Concentrations</b>	0.5 mg/ml		
<b>Immunogen</b>	Human IgA.		

**External Database****Links****UniProt:**[P01876](#)    [Related reagents](#)[P01877](#)    [Related reagents](#)**Entrez Gene:**[3493](#) IGHA1    [Related reagents](#)[3494](#) IGHA2    [Related reagents](#)**RRID**

AB\_616850

**Specificity**

**Goat anti Human IgA antibody** recognizes the heavy chain of human IgA as demonstrated by ELISA and Flow cytometry. IgA is the major category of antibody present in mucosal secretions and acts as a first line of defense against inhaled or ingested pathogens at the mucosal interface. IgA is also present in serum where it acts to eliminate pathogens via interaction with FcαR, also known as CD89. ([Woof et al. 2004](#)). Biotinylated Goat anti Human IgA antibody (**205008**) has been cross-adsorbed against human IgD, IgG and IgM to ensure class specificity.

**References**

1. Wright, A.K. et al. (2012) Human Nasal Challenge with *Streptococcus pneumoniae* Is Immunising in the Absence of Carriage. [PLoS Pathog. 8: e1002622](#).
2. Pennington, S.H. et al. (2016) Polysaccharide-specific Memory B-cells Predict Protection Against Experimental Human Pneumococcal Carriage. [Am J Respir Crit Care Med. Jul 12. \[Epub ahead of print\]](#)
3. Turner, P.J. et al. (2020) Differences in nasal immunoglobulin A responses to influenza vaccine strains after live attenuated influenza vaccine (LAIV) immunization in children. [Clin Exp Immunol. 199 \(2\): 109-118](#).

**Storage**

+4°C DO NOT FREEZE

This product is photosensitive and should be protected from light.

**Guarantee**

Guaranteed until date of expiry. Please see product label.

**Health And Safety Information**

Material Safety Datasheet documentation #10331 available at:  
<https://www.bio-rad-antibodies.com/SDS/205009>  
10331

**Regulatory**

For research purposes only

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)[HUMAN SEROBLOCK \(BUF070B\)](#)**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M386895:210601'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)