

Datasheet: 2083-4955

BATCH NUMBER 165354

Description:	MOUSE ANTI CHONDROITIN SULPHATE
Specificity:	CHONDROITIN SULPHATE
Format:	Ascites
Product Type:	Monoclonal Antibody
Clone:	CS-56
Isotype:	IgM
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/100 - 1/400
Immunofluorescence	▪			1/100 - 1/400

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species

Chicken

Species Cross Reactivity

Reacts with: Bovine, Human, Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Ascitic Fluid - raw

Preservative Stabilisers

0.1% Sodium Azide (NaN₃)

Immunogen

Ventral membranes of chicken gizzard fibroblasts.

RRID

AB_620552

Specificity	<p>Mouse anti Chondroitin sulphate antibody, clone CS-56 recognizes chondroitin sulphate types A and C (not B, dermatan sulphate).</p> <p>Mouse anti Chondroitin sulphate antibody, clone CS-56 gives positive staining in bovine mammary gland epithelial cells, and also in human, chicken, and mouse fibroblast or tissues.</p>
References	<p>1. Avnur, Z. <i>et al</i> (1985) Spatial interrelationships between proteoglycans and extracellular matrix proteins in cell cultures. Exp Cell Res. 158:321-32.</p> <p>2. Avnur, Z.<i>et al</i> (1984) Immunocytochemical Localization of Native Chondroitin-Sulphate in Tissues and Cultured Cells Using Specific Monoclonal Antibody. Cell 38: 811-822.</p>
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10342 available at: https://www.bio-rad-antibodies.com/SDS/2083-4955</p> <p>10342</p>
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M381822:210512'

Printed on 29 Aug 2024