

Datasheet: 2083-4955

**BATCH NUMBER 180221**

<b>Description:</b>	MOUSE ANTI CHONDROITIN SULPHATE
<b>Specificity:</b>	CHONDROITIN SULPHATE
<b>Format:</b>	Ascites
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	CS-56
<b>Isotype:</b>	IgM
<b>Quantity:</b>	0.1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/100 - 1/400
Immunofluorescence	▪			1/100 - 1/400

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

### Target Species

Chicken

### Species Cross Reactivity

Reacts with: Bovine, Human, Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Ascitic Fluid - raw

### Preservative Stabilisers

0.1% Sodium Azide (NaN<sub>3</sub>)

### Immunogen

Ventral membranes of chicken gizzard fibroblasts.

### RRID

AB\_620552

**Specificity** **Mouse anti Chondroitin sulphate antibody, clone CS-56** recognizes chondroitin sulphate types A and C (not B, dermatan sulphate).

Mouse anti Chondroitin sulphate antibody, clone CS-56 gives positive staining in bovine mammary gland epithelial cells, and also in human, chicken, and mouse fibroblast or tissues.

---

**References**

1. Avnur, Z. *et al* (1985) Spatial interrelationships between proteoglycans and extracellular matrix proteins in cell cultures. [Exp Cell Res. 158:321-32.](#)
2. Avnur, Z.*et al* (1984) Immunocytochemical Localization of Native Chondroitin-Sulphate in Tissues and Cultured Cells Using Specific Monoclonal Antibody. [Cell 38: 811-822.](#)

---

**Storage**

Store at +4°C or at -20°C if preferred.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10342 available at: <https://www.bio-rad-antibodies.com/SDS/2083-4955>  
10342

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M362971:200528'

**Printed on 15 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)