

Datasheet: 2150-0020

BATCH NUMBER 157443

Description:	RABBIT ANTI HUMAN COLLAGEN I
Specificity:	COLLAGEN I
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/10 - 1/50
Immunohistology - Paraffin (1)	▪			1/10 - 1/50
ELISA	▪			1/500 - 1/3000
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

(1) This product requires antigen retrieval using 1mg/ml pepsin in 0.5M acetic acid for 2 hours at 37°C and multiple buffer washes prior to staining of paraffin sections.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	<p>Reconstitute with 0.5ml distilled water</p> <p>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 1% BSA and 0.09% sodium azide is recommended.</p>
Preparation	Purified IgG prepared by affinity chromatography on immobilized Collagen type I
Buffer Solution	<p>0.15M Sodium Chloride</p> <p>0.01M Sodium Phosphate</p>

Preservative	1% Mannitol
Stabilisers	1% Dextran
Immunogen	Collagen type I from human placenta.
External Database Links	<p>UniProt:</p> <p>P02452 Related reagents</p> <p>P08123 Related reagents</p> <p>Entrez Gene:</p> <p>1277 COL1A1 Related reagents</p> <p>1278 COL1A2 Related reagents</p>
RRID	AB_620307
Specificity	Rabbit anti Human collagen I antibody recognizes human Collagen I, reacting with both native and heat denatured proteins, and has been cross-absorbed against Collagen types III, IV and V.
References	<ol style="list-style-type: none"> Evans, M.J. <i>et al.</i> (2002) Fibroblast growth factor-2 during postnatal development of the tracheal basement membrane zone. Am J Physiol Lung Cell Mol Physiol. 283 (6): L1263-70. Kim, H.J. & Im, G.I. (2010) The effects of ERK1/2 inhibitor on the chondrogenesis of bone marrow- and adipose tissue-derived multipotent mesenchymal stromal cells. Tissue Eng Part A. 16: 851-60. Lee, J.S. and Im, G.I. (2010) Influence of chondrocytes on the chondrogenic differentiation of adipose stem cells. Tissue Eng Part A. 16: 3569-77. Kim, Y.J. <i>et al.</i> (2008) PTHrP promotes chondrogenesis and suppresses hypertrophy from both bone marrow-derived and adipose tissue-derived MSCs. Biochem Biophys Res Commun. 373: 104-8 Li M <i>et al.</i> (2013) Short periods of cyclic mechanical strain enhance triple-supplement directed osteogenesis and bone nodule formation by human embryonic stem cells <i>in vitro</i>. Tissue Eng Part A. 19 (19-20): 2130-7. Hoffman, J.K. <i>et al.</i> (2015) Articular cartilage repair using marrow stimulation augmented with a viable chondral allograft: 9-month postoperative histological evaluation. Case Rep Orthop. 2015: 617365.
Storage	<p>Prior to reconstitution store at +4°C.</p> <p>After reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	<p>Material Safety Datasheet documentation #20472 available at: https://www.bio-rad-antibodies.com/SDS/2150-0020</p> <p>20472</p>

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M374982:201210'

Printed on 27 Sep 2024