

Datasheet: 2150-0020

BATCH NUMBER 160084

Description:	RABBIT ANTI HUMAN COLLAGEN I
Specificity:	COLLAGEN I
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/10 - 1/50
Immunohistology - Paraffin (1)	▪			1/10 - 1/50
ELISA	▪			1/500 - 1/3000
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

(1)This product requires antigen retrieval using 1mg/ml pepsin in 0.5M acetic acid for 2 hours at 37°C and multiple buffer washes prior to staining of paraffin sections.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.5ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 1% BSA and 0.09% sodium azide is recommended.
Preparation	Purified IgG prepared by affinity chromatography on immobilized Collagen type I
Buffer Solution	0.1M Sodium Chloride 0.01M Sodium Phosphate

0.01M Sodium Borate

Preservative 1% Mannitol
Stabilisers 1% Dextran

Immunogen Collagen type I from human placenta.

External Database Links

UniProt:

[P02452](#) [Related reagents](#)

[P08123](#) [Related reagents](#)

Entrez Gene:

[1277](#) COL1A1 [Related reagents](#)

[1278](#) COL1A2 [Related reagents](#)

RRID AB_620307

Specificity **Rabbit anti Human collagen I antibody** recognizes human Collagen I, reacting with both native and heat denatured proteins, and has been cross-absorbed against Collagen types III, IV and V.

References

1. Evans, M.J. *et al.* (2002) Fibroblast growth factor-2 during postnatal development of the tracheal basement membrane zone. [Am J Physiol Lung Cell Mol Physiol. 283 \(6\): L1263-70.](#)
 2. Kim, Y.J. *et al.* (2008) PTHrP promotes chondrogenesis and suppresses hypertrophy from both bone marrow-derived and adipose tissue-derived MSCs. [Biochem Biophys Res Commun. 373: 104-8](#)
 3. Kim, H.J. & Im, G.I. (2010) The effects of ERK1/2 inhibitor on the chondrogenesis of bone marrow- and adipose tissue-derived multipotent mesenchymal stromal cells. [Tissue Eng Part A. 16: 851-60.](#)
 4. Lee, J.S. and Im, G.I. (2010) Influence of chondrocytes on the chondrogenic differentiation of adipose stem cells. [Tissue Eng Part A. 16: 3569-77.](#)
 5. Li M *et al.* (2013) Short periods of cyclic mechanical strain enhance triple-supplement directed osteogenesis and bone nodule formation by human embryonic stem cells *in vitro*. [Tissue Eng Part A. 19 \(19-20\): 2130-7.](#)
 6. Hoffman, J.K. *et al.* (2015) Articular cartilage repair using marrow stimulation augmented with a viable chondral allograft: 9-month postoperative histological evaluation. [Case Rep Orthop. 2015: 617365.](#)
 7. Segreto, F. *et al.* (2020) The use of acellular porcine dermis, hyaluronic acid and polynucleotides in the treatment of cutaneous ulcers: Single blind randomised clinical trial. [Int Wound J. 17 \(6\): 1702-8.](#)
-

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #20482 available at:
<https://www.bio-rad-antibodies.com/SDS/2150-0020>
20482

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

 Tel: +44 (0)1865 852 700

 Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

 Tel: +49 (0) 89 8090 95 21

 Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M391125:211008'

Printed on 27 Sep 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)