

Datasheet: 2150-0030

**BATCH NUMBER 157150**

<b>Description:</b>	NATIVE HUMAN COLLAGEN I
<b>Name:</b>	COLLAGEN I
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified protein from human placenta - lyophilized
<b>Reconstitution</b>	Use 0.05 - 0.5 M acetic acid, pH 2.5 at 4°C.
<b>Preparation</b>	Collagen from washed dissected tissue was dissolved in acetic acid after mild pepsin treatment. Collagen type I was then purified using differential salt precipitation.
<b>Preservative Stabilisers</b>	None present

### External Database Links

#### UniProt:

[P02452](#)    [Related reagents](#)

[P08123](#)    [Related reagents](#)

#### Entrez Gene:

[1277](#) COL1A1    [Related reagents](#)

[1278](#) COL1A2    [Related reagents](#)

---

**Product Information** Collagen is located in the extracellular matrix of connective tissues. It is part of the interacting network of proteoglycans and proteins that provides a structural framework for both soft and calcified connective tissues. Type I collagen (95 kDa) is found in bone, cornea, skin and tendon. Mutations in the encoding gene are associated with osteogenesis imperfecta, Ehlers Danlos syndrome, and idiopathic osteoporosis.

---

**Protein Molecular Weight** 300 kDa

---

**Purity** >90%

---

**Storage** Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

---

**Guarantee** Guaranteed until date of expiry. Please see product label.

---

**Acknowledgements** The human source was tested and found negative for HBsAg and HIV-1/2 antibodies.

---

**Health And Safety Information** Material Safety Datasheet documentation #10268 available at:  
<https://www.bio-rad-antibodies.com/SDS/2150-0030>  
10268

Donor material tested and found negative for HBsAg, HCV, and HIV1 and 2 antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M359497:190924'

**Printed on 25 Mar 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)