

Datasheet: 2150-0515

BATCH NUMBER 165247

Description:	NATIVE BOVINE COLLAGEN I
Name:	COLLAGEN I
Format:	Purified
Product Type:	Purified Protein
Quantity:	10 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified protein from calf tendon - lyophilized
Reconstitution	Use 0.5M acetic acid, pH 3.0. Dissolved collagen have to be thoroughly homogenized in a Potter after addition of the solute and left overnight in the coldroom to dissolve. Dissolved collagen retains immunologic properties of native collagen type I. Thermal denaturation (5 minutes at 100°C) converts dissolved collagen to gelatin.
Preparation	Bovine collagen type I was isolated from calf tendons. Isolation steps included extraction into dilute acetic acid after pepsin treatment, followed by multiple steps of differential salt precipitation and ion-exchange chromatography.
Buffer Solution	Essentially salt free
Preservative Stabilisers	None present
External Database Links	UniProt:

[P02465](#) [Related reagents](#)

[P02453](#) [Related reagents](#)

Entrez Gene:

[282188](#) COL1A2 [Related reagents](#)

[282187](#) COL1A1 [Related reagents](#)

Product Information Collagen type I 90%
Collagen type III 10%
Non-collagenous proteins <0.5%

Protein Molecular Weight ~300 kDa

Storage Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Avoid repeated freezing and thawing as this may denature the antibody.
Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10268 available at:
<https://www.bio-rad-antibodies.com/SDS/2150-0515>
10268

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350171:190307'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)