

Datasheet: 2150-2203

**BATCH NUMBER 165646**

<b>Description:</b>	MOUSE ANTI HUMAN COLLAGEN I/II/III
<b>Specificity:</b>	COLLAGEN I/II/III
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	NLI/22
<b>Isotype:</b>	IgM
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/5 - 1/20
ELISA	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgM - lyophilized
<b>Reconstitution</b>	Use 1 ml sterile distilled water to give a concentration of 0.1 mg/ml.
<b>Preparation</b>	Gel Filtration chromatography from ascites
<b>Buffer Solution</b>	0.1M Sodium Chloride 0.01M Sodium Phosphate 0.01M Sodium Borate
<b>Preservative Stabilisers</b>	1% Mannitol 1% Dextran
<b>Immunogen</b>	Human placental collagen type 1.

**External Database  
Links**

**UniProt:**

<a href="#">P02452</a>	<a href="#">Related reagents</a>
<a href="#">P08123</a>	<a href="#">Related reagents</a>
<a href="#">P02458</a>	<a href="#">Related reagents</a>
<a href="#">P02461</a>	<a href="#">Related reagents</a>

**Entrez Gene:**

<a href="#">1277</a>	COL1A1	<a href="#">Related reagents</a>
<a href="#">1278</a>	COL1A2	<a href="#">Related reagents</a>
<a href="#">1280</a>	COL2A1	<a href="#">Related reagents</a>
<a href="#">1281</a>	COL3A1	<a href="#">Related reagents</a>

---

**RRID**

AB\_617237

---

**Specificity**

**Rabbit anti Human collagen I/II/III antibody** reacts with collagens type I, II and III. It does not react with human collagen types IV, V or VI, human plasma proteins or albumin from various species.

---

**Storage**

Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

Guaranteed until date of expiry. Please see product label.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #20482 available at:  
<https://www.bio-rad-antibodies.com/SDS/2150-2203>  
20482

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgM (STAR138...) [Alk. Phos.](#)

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)

**North & South  
America**

Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M405726:220916'

**Printed on 20 Mar 2024**