

Datasheet: 2150-2206

BATCH NUMBER 170033

Description:	RABBIT ANTI HUMAN COLLAGEN I/II/III/IV/V
Specificity:	COLLAGEN I/II/III/IV/V
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/10 - 1/40
Immunohistology - Paraffin			▪	
ELISA	▪			1/2000 - 1/8000
Western Blotting			▪	
Immunofluorescence	▪			1/10 - 1/40

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.5ml distilled water.
Preparation	Purified Ig prepared by affinity chromatography on immobilized collagen types I-V
Buffer Solution	Borate buffered saline
Preservative	1% Mannitol
Stabilisers	1% Dextran
Approx. Protein Concentrations	IgG concentration 0.1mg/ml after reconstitution.

Immunogen

A mixture of human collagen I, II, III, IV and V from human placenta.

**External Database
Links****UniProt:**

[P02452](#) [Related reagents](#)
[P08123](#) [Related reagents](#)
[P02458](#) [Related reagents](#)
[P02461](#) [Related reagents](#)
[P02462](#) [Related reagents](#)
[P08572](#) [Related reagents](#)
[Q01955](#) [Related reagents](#)
[P53420](#) [Related reagents](#)
[P29400](#) [Related reagents](#)
[Q14031](#) [Related reagents](#)
[P20908](#) [Related reagents](#)
[P05997](#) [Related reagents](#)
[P25940](#) [Related reagents](#)

Entrez Gene:

[1277](#) COL1A1 [Related reagents](#)
[1278](#) COL1A2 [Related reagents](#)
[1280](#) COL2A1 [Related reagents](#)
[1281](#) COL3A1 [Related reagents](#)
[1282](#) COL4A1 [Related reagents](#)
[1284](#) COL4A2 [Related reagents](#)
[1285](#) COL4A3 [Related reagents](#)
[1286](#) COL4A4 [Related reagents](#)
[1287](#) COL4A5 [Related reagents](#)
[1288](#) COL4A6 [Related reagents](#)
[1289](#) COL5A1 [Related reagents](#)
[1290](#) COL5A2 [Related reagents](#)
[50509](#) COL5A3 [Related reagents](#)

RRID

AB_620306

Specificity

Rabbit anti Human collagen I/II/III/IV/V antibody recognizes multiple collagen types and exhibits the following reactivities in ELISA with both native and denatured human collagen types:

I & III 100%
IV & V 90%
II 80%

Human fibronectin, albumin and immunoglobulins all <0.5%

Storage This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #20482 available at: <https://www.bio-rad-antibodies.com/SDS/2150-2206>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M429172:240320'

Printed on 04 Mar 2026