

Datasheet: 2150-2555

**BATCH NUMBER 157444**

<b>Description:</b>	RABBIT ANTI PIG COLLAGEN I/III
<b>Specificity:</b>	COLLAGEN I/III
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/10 - 1/40
ELISA	▪			1/500 - 1/2000
Western Blotting			▪	
Immunofluorescence	▪			1/10 - 1/40

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute the vial (50 µg) with 500 µl of sterile distilled water or saline.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on porcine collagen types I & III
<b>Buffer Solution</b>	0.15M Sodium Chloride 0.01M Sodium Phosphate
<b>Preservative Stabilisers</b>	1% Mannitol 1% Dextran
<b>Immunogen</b>	Porcine collagen, types 1 and 3 from skin.

RRID AB\_618709

---

**Specificity** **Rabbit anti Porcine collagen I/III antibody** reacts with native and heat denatured porcine collagen type 1 and 3.  
Less than 0.5% reactivity is observed with porcine albumin and immunoglobulins.

---

**Storage** Prior to reconstitution store at +4°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.

---

**Guarantee** Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information** Material Safety Datasheet documentation #20472 available at:  
<https://www.bio-rad-antibodies.com/SDS/2150-2555>  
20472

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M374981:201210'

Printed on 12 Aug 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)