

Datasheet: 2150-2555

**BATCH NUMBER 159095**

<b>Description:</b>	RABBIT ANTI PIG COLLAGEN I/III
<b>Specificity:</b>	COLLAGEN I/III
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.5 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			1/10 - 1/40
ELISA	■			1/500 - 1/2000
Western Blotting			■	
Immunofluorescence	■			1/10 - 1/40

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute the vial (50 µg) with 500 µl of sterile distilled water or saline.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on porcine collagen types I & III
<b>Buffer Solution</b>	0.1M Sodium Chloride 0.01M Sodium Phosphate 0.01M Sodium Borate
<b>Preservative Stabilisers</b>	1% Mannitol 1% Dextran
<b>Immunogen</b>	Porcine collagen, types 1 and 3 from skin.

<b>RRID</b>	AB_618709
<b>Specificity</b>	<b>Rabbit anti Porcine collagen I/III antibody</b> reacts with native and heat denatured porcine collagen type 1 and 3. Less than 0.5% reactivity is observed with porcine albumin and immunoglobulins.
<b>Storage</b>	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20482 available at: <a href="https://www.bio-rad-antibodies.com/SDS/2150-2555">https://www.bio-rad-antibodies.com/SDS/2150-2555</a> 20482
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M380410:210511'

Printed on 20 Mar 2024