

Datasheet: 2330-4809

**BATCH NUMBER 164443**

<b>Description:</b>	MOUSE ANTI CORTISOL
<b>Specificity:</b>	CORTISOL
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	C13 (BGN/C13)
<b>Isotype:</b>	IgG3
<b>Quantity:</b>	0.2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/50 - 1/1500
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

### Target Species

Broad

### Species Cross Reactivity

Reacts with: Human  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

### Buffer Solution

Phosphate buffered saline

<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )										
<b>Carrier Free</b>	Yes										
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml										
<b>Immunogen</b>	Cortisol conjugated to bovine serum albumin.										
<b>RRID</b>	AB_618933										
<b>Specificity</b>	<p><b>Mouse anti cortisol monoclonal antibody, clone C13</b> recognizes cortisol, a glucocorticoid hormone synthesized from cholesterol by the adrenal cortex. It is released in response to pituitary adenocorticotrophic hormone (ACTH) which is, in turn, regulated by corticotrophin releasing factor (CRF). High levels of cortisol act to regulate the release of ACTH and CRF via negative feedback inhibition.</p> <p>The release of cortisol stimulates the conversion of proteins to carbohydrates, raises blood sugar levels and promotes glycogen storage in the liver.</p> <p>The following cross reactivities were observed:</p> <table border="0"> <tr> <td>Cortisol</td> <td>100%</td> </tr> <tr> <td>Prednisolone</td> <td>30%</td> </tr> <tr> <td>Cortisone</td> <td>25%</td> </tr> <tr> <td>Corticosterone</td> <td>10%</td> </tr> <tr> <td>11-desoxycortisol</td> <td>5%</td> </tr> </table> <p>No reactivity was observed with 17-a-OH-progesterone.</p>	Cortisol	100%	Prednisolone	30%	Cortisone	25%	Corticosterone	10%	11-desoxycortisol	5%
Cortisol	100%										
Prednisolone	30%										
Cortisone	25%										
Corticosterone	10%										
11-desoxycortisol	5%										
<b>References</b>	<ol style="list-style-type: none"> <li>1. Arya, S.K. <i>et al.</i> (2010) Antibody functionalized interdigitated micro-electrode (IDmicroE) based impedimetric cortisol biosensor. <a href="#">Analyst. 135 (8): 1941-6.</a></li> <li>2. Arya, S.K. <i>et al.</i> (2010) Dithiobis(succinimidyl propionate) modified gold microarray electrode based electrochemical immunosensor for ultrasensitive detection of cortisol. <a href="#">Biosens Bioelectron. 25: 2296-301.</a></li> <li>3. Arya, S.K. <i>et al.</i> (2011) Polyaniline protected gold nanoparticles based mediator and label free electrochemical cortisol biosensor. <a href="#">Biosens Bioelectron. 28 (1): 166-73.</a></li> <li>4. Kaushik, A. <i>et al.</i> (2013) Mediator and label free estimation of stress biomarker using electrophoretically deposited Ag@AgO-polyaniline hybrid nanocomposite. <a href="#">Biosens Bioelectron. 50: 35-41.</a></li> </ol>										
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>										
<b>Guarantee</b>	12 months from date of despatch										

---

**Health And Safety Information**      Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/2330-4809>  
10040

---

**Regulatory**                      For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)                      [HRP](#)  
Rabbit Anti Mouse IgG (STAR12...)                      [RPE](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)  
Goat Anti Mouse IgG (STAR76...)                      [RPE](#)  
Rabbit Anti Mouse IgG (STAR13...)                      [HRP](#)  
Goat Anti Mouse IgG (STAR70...)                      [FITC](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),  
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),  
[FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...)                      [FITC](#)  
Goat Anti Mouse IgG (Fc) (STAR120...)                      [FITC](#), [HRP](#)

**North & South America**      Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M405264:220916'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)