

Datasheet: 302009

Description:	GOAT ANTI RAT IgM:RPE
Specificity:	IgM
Format:	RPE
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

RRID AB_616872

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat - 1/5

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Rat

Product Form Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578

Antiserum Preparation Antiserum to rat IgM was raised by repeated immunisation of goats with highly purified antigen. The preparation was purified by affinity chromatography on rat IgM covalently linked to agarose.

Buffer Solution Phosphate buffered saline

Preservative 0.1% Sodium azide
Stabilisers Stabilising agent

Approx. Protein Concentrations IgG concentration 0.5 mg/ml

Immunogen Normal rat IgM.

Specificity **Goat anti Rat IgM antibody** recognizes rat IgM as demonstrated by ELISA. This antibody has been cross-absorbed against rat IgG, but may react with IgM from other species.

Flow Cytometry Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.

References 1. Brissette, C.A. *et al.* (2012) Acetate supplementation reduces microglia activation and brain interleukin-1 β levels in a rat model of Lyme neuroborreliosis. [J Neuroinflammation. 9: 249.](#)

Storage +4°C DO NOT FREEZE
This product is photosensitive and should be protected from light.
This product should be stored undiluted.
Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10331 available at:
10331: <https://www.bio-rad-antibodies.com/uploads/MSDS/10331.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M350184:190307'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)